

NRC Publications Archive Archives des publications du CNRC

Tessier sequential extraction on 17 elements from three marine sediment certified reference materials (HISS-1, MESS-4, and PACS-3)
Kumkrong, Paramee; Mihai, Ovi; Mercier, Patrick H. J.; Pihilligawa Gedara, Indu; Tyo, Daniel D.; Mester, Zoltán

This publication could be one of several versions: author's original, accepted manuscript or the publisher's version. / La version de cette publication peut être l'une des suivantes : la version prépublication de l'auteur, la version acceptée du manuscrit ou la version de l'éditeur.

For the publisher's version, please access the DOI link below. / Pour consulter la version de l'éditeur, utilisez le lien DOI ci-dessous.

Publisher's version / Version de l'éditeur:

<https://doi.org/10.1007/s00216-020-03063-z>

Analytical and Bioanalytical Chemistry, 2020-11-24

NRC Publications Archive Record / Notice des Archives des publications du CNRC :

<https://nrc-publications.canada.ca/eng/view/object/?id=e0281ead-bdd3-4b22-a2cb-38344aa397a1>

<https://publications-cnrc.canada.ca/fra/voir/objet/?id=e0281ead-bdd3-4b22-a2cb-38344aa397a1>

Access and use of this website and the material on it are subject to the Terms and Conditions set forth at

<https://nrc-publications.canada.ca/eng/copyright>

READ THESE TERMS AND CONDITIONS CAREFULLY BEFORE USING THIS WEBSITE.

L'accès à ce site Web et l'utilisation de son contenu sont assujettis aux conditions présentées dans le site

<https://publications-cnrc.canada.ca/fra/droits>

LISEZ CES CONDITIONS ATTENTIVEMENT AVANT D'UTILISER CE SITE WEB.

Questions? Contact the NRC Publications Archive team at

PublicationsArchive-ArchivesPublications@nrc-cnrc.gc.ca. If you wish to email the authors directly, please see the first page of the publication for their contact information.

Vous avez des questions? Nous pouvons vous aider. Pour communiquer directement avec un auteur, consultez la première page de la revue dans laquelle son article a été publié afin de trouver ses coordonnées. Si vous n'arrivez pas à les repérer, communiquez avec nous à PublicationsArchive-ArchivesPublications@nrc-cnrc.gc.ca.

Tessier sequential extraction on 17 elements from three marine sediment certified reference materials (HISS-1, MESS-4 and PACS-3)

Paramee Kumkrong*, Ovi Mihai, Patrick H.J. Mercier, Indu Gedara Pihilligawa, Daniel D. Tyo, and Zoltan Mester

National Research Council Canada, 1200 Montreal Road Ottawa, Ontario Canada, K1A 0R6

*Corresponding author : paramee.kumkrong@nrc-crnc.gc.ca

ORCID : 0000-0002-9281-9974

Abstract

The four-step Tessier sequential extraction procedure is a well-known approach used for environmental and geochemical studies in soil and sediments. However, a lack of reference materials limits its use making implementation and quality control cumbersome. This study applied Tessier sequential extraction to three globally used marine sediment certified reference materials (CRMs) including HISS-1, MESS-4 and PACS-3 with varying levels of contamination.

The study analyzed the distribution of 17 elements throughout the extraction phases. Overall, the percent recovery (sum of steps vs total metal concentration) of all analyzed elements in Tessier extraction was $92 \% \pm 40 \%$ in HISS-1, $101 \% \pm 12 \%$ in MESS-4 and $102 \% \pm 10 \%$ in PACS-3. The observed uncertainty of the individual elemental concentrations averaged of 13 %, which compares favorably with the 16 % target uncertainty derived from Horwitz equation. The reference data set produced here using the Tessier sequential extraction procedure will serve as a quality control and method development tool for laboratories.

Keywords: Distribution. Certified reference material. Sequential extraction. Recovery. Uncertainty.

24	Table of contents	
25	Abstract	1
26	Abbreviation	3
27	1 Introduction.....	5
28	2. Materials and methods	6
29	2.1 Method for sediment and residue digestion	7
30	2.2 Tessier sequential extraction procedure.....	8
31	2.3 Data manipulation	10
32	3 Results.....	11
33	3.1 Total concentrations of elements in HISS-1, MESS-4, and PACS-3	11
34	3.2 Concentrations of elements by Tessier sequential extraction	11
35	4 Discussion	20
36	4.1 Total concentrations of elements in HISS-1, MESS-4, and PACS-3	20
37	4.2 Distribution profile of elements from Tessier sequential extraction	21
38	4.2.1 HISS-1	21
39	4.2.2 MESS-4.....	21
40	4.2.3 PACS-3	22
41	4.3 Analytical performance of Tessier sequential extraction on three marine sediment CRMs.....	22
42	4.3.1 Assessment of precision.....	22
43	4.3.2 Assessment of accuracy	25
44	5 Conclusions.....	27
45	Acknowledgement	27
46	Compliance with ethical standards	28
47	Conflict of interest	28
48	References.....	28

1
2
3 51 **Abbreviation**
4
5

6		
7		
8	°C	Degree Celsius
9		
10	µm	Micrometer
11		
12	BCR	Community Bureau of Reference
13		
14	CI	Confidence interval
15		
16	CRC	Collision reaction cell
17		
18	CRM(s)	Certified reference material(s)
19		
20	DIW	Deionized water
21		
22	eV	Electron volt
23		
24	Ext sol#	Extraction solution
25		
26	g	Gram
27		
28	h	Hour
29		
30	ICP	Inductively coupled plasma
31		
32	ISO	International Organization for Standardization
33		
34	L	Litre
35		
36	log	Logarithm
37		
38	mol L ⁻¹	Mole per litre
39		
40	mg	Milligram
41		
42	min	Minute
43		
44	mL	Millilitre
45		
46	MS	Mass spectrometry
47		
48	MΩ-cm	Megohm-centimeter
49		
50	na	Not applicable
51		
52	nd	Not detected
53		
54	nr	Not reported
55		
56	NRC	National Research Council Canada
57		

OES	Optical emission spectrometry
rpm	Revolutions per minute
RSD	Relative standard deviation
SD	Standard deviation
S_R	Reproducibility standard deviation
u_{std}	Standard uncertainty
U	Expanded uncertainty
v	Volume
w	Weight

52

53

1 Introduction

The association between metals bound to sediments is of interest because their potential release could pose significant environmental challenges (1,2). Marine sediments are mainly composed of biogenous, hydrogenous, and lithogenous sediments (3), and also receive foreign materials from inland activities such as weathering and anthropogenic inputs. A high accumulation of toxic metals in marine sediments is widely observed in areas of high human activity, for example harbors (4–8).

In general, assessment of sediment quality is achieved by quantitative analysis of total toxic metal content (9–11). However, not all measured total metal content in sediments is actually posing an environmental hazard, so better approaches that estimate the labile metal fraction and the conditions of their release are needed. The importance of such investigations brought the development of methodologies based on fractionating toxic metals from solid phases by partial dissolution / desorption / extraction (12–17).

Among sequential extraction procedures, the Tessier and the Community Bureau of Reference (BCR) sequential extractions are most commonly employed (18,19). The publications citing Tessier and BCR sequential extraction in Scopus database (<https://www.scopus.com>, accessed on June 7, 2020) showed results of 3,625 and 2,332, respectively. Tessier expressed significant growth until it reached a plateau in 2011, however, BCR extraction method has continued gaining popularity (**Fig.1**).

Tessier sequential extraction procedure was published in 1979 for the determination of Cd, Cr, Cu, Fe, Mn, Ni, Pb, and Zn in river sediments from Quebec, Canada (18). The BCR method was developed with the objective to harmonize the sequential extraction procedure and be used in conjunction with the production of BCR 601 and BCR 701 CRMs (river sediment reference materials) that certified Cd, Cr, Cu, Ni, Pb, and Zn (19,20). There is a plethora of data available from the BCR on method characteristics (21–26) but the availability of such information is much limited for the Tessier extraction approach using reference materials.




Figure 1 : The number of publications of Tessier and BCR sequential extraction procedure since 1979 to 2020 (Source: Scopus database, searched under keywords: “sequential extraction BCR” and “sequential extraction Tessier” (<https://www.scopus.com>, accessed on June 07, 2020))

Considering the lack of CRMs and reference data on precision and accuracy upon application of a given extraction method, we have performed the Tessier extraction procedure on 17 elements (As, Be, Cd, Cr, Cu, Fe, Li, Mn, Ni, Pb, Sb, Se, Sn, Sr, U, V, and Zn) in three marine sediment CRMs produced by National Research Council Canada (NRC). The three CRMs included: HISS-1, a clean sediment; MESS-4, a moderately-contaminated sediment; and PACS-3, a highly-contaminated marine sediment (27–29). These were certified for total concentrations of numerous elements and have been used as quality control and for method validation purpose in analytical laboratories worldwide.

2. Materials and methods

HISS-1 CRM was sampled from Hibernia shelf, off the coast of Newfoundland, Canada (27), MESS-4 from Beaufort Sea, Canadian Arctic (28) and PACS-3 from Esquimalt harbor, British Columbia, Canada (29). The CRMs were sieved to particle sizes of less than 125 μm . Homogeneity and stability was assessed. HISS-1 CRM has served as marine sediment reference materials since 1997, PACS-3 since 2013 and MESS-4 since 2014. The three marine sediment CRMs were used “as is” and the moisture content was corrected to obtain concentration based on dry weight.

Laboratory ware, chemicals, and instruments are listed in Table S1 of Electronic Supplementary Material. The extraction solution (Ext sol #) and its preparation in each step are provided in **Table 1**.

Table 1: Extraction solutions for Tessier sequential extraction

Ext sol #	Reagent	Reagent preparation
1	1 mol L ⁻¹ sodium acetate solution (NaOAc), pH 8.2	Dissolve 41 g of NaOAc in 450 mL of deionized water (DIW), adjust to pH 8.2 using acetic acid (AcOH) and dilute to 500 mL with DIW
2	1 mol L ⁻¹ NaOAc, pH 5.0	Same as Ext sol# 1 but adjust to pH 5 using AcOH
3	0.04 mol L ⁻¹ hydroxylammonium chloride solution (NH ₂ OH.HCl) in 25 % (v/v) AcOH	Dissolve 1.4 g of NH ₂ OH.HCl in 500 mL of 25 % (v/v) AcOH
4	0.02 mol L ⁻¹ nitric acid (HNO ₃)	- Add 6.3 mL of 69 % (w/w) HNO ₃ and make up to 1,000 mL with DIW to obtain 0.1 mol L ⁻¹ HNO ₃ - Pipette 20 mL of 0.1 mol L ⁻¹ HNO ₃ and make up to 100 mL with DIW to obtain 0.02 mol L ⁻¹ HNO ₃
	30 % (w/v) hydrogen peroxide (H ₂ O ₂) at pH 2	Add few drops of 69% (w/w) HNO ₃ to 30% (w/v) H ₂ O ₂ and check pH with pH strip to obtain pH 2
	3.2 mol L ⁻¹ ammonium acetate solution (NH ₄ OAc) in 20 % (v/v) HNO ₃	Dissolve 24.7 g of NH ₄ OAc in 100 mL of 20 % (v/v) of HNO ₃

2.1 Method for sediment and residue digestion

The following digestion method is used for total sediment digestion of the raw material CRMs and reported results in Section 3.1 Total concentrations of elements and Section 3.2 Concentrations of elements by Tessier sequential extraction in Table 2 to 4 for the residues. Two laboratories (Lab A and Lab B) used their routine method for digestion as described by Kumkrong et.al. (21). In brief Lab A used 0.25 g of sample and a mixture of HNO₃, hydrochloric acid (HCl) and hydrogen fluoride (HF). The excess amount of HF was destroyed by adding

1
2
3 101 boric acid (H_3BO_3). After the microwave digestion, the samples were placed in a hot block, evaporated to obtain 5 mL
4
5 102 volume and diluted to 50 g using 2 % (v/v) HNO_3 .
6

7
8 103 Lab B used 0.25 g of raw material or 1 g of residue and digested using a mixture of HNO_3 , HF, and H_2O_2 . The opened
9
10 104 digestion was performed after microwave digestion using a mixture of HNO_3 and HCl on a hot block and evaporated
11
12 105 to approximately 1 mL and the digestate was diluted to 100 g using DIW.
13

14 106 The elemental analysis of digestates was conducted using inductively coupled plasma-optical emission
15
16 107 spectrometry (ICP-OES) and inductively coupled plasma-mass spectrometry (ICP-MS). The conditions of
17
18 108 ICP-OES and ICP-MS are stated in the Electronic Supplementary Material Table S2.
19

20 21 109 **2.2 Tessier sequential extraction procedure**

22
23 110 The extraction procedure is presented in **Fig 2.** and detailed in Table S3 of Electronic Supplementary Material. The
24
25 111 extraction tube was weighed before and after adding a new extraction solution for a weight correction. The extraction
26
27 112 steps were labeled as Tessier step 1 - exchangeable (T # 1), step 2 - bound to carbonates (T # 2), step 3 - bound to Fe-
28
29 113 Mn oxides (T # 3), step 4 - bound to organic matter and sulfide minerals (T # 4), and residue. The fractions were
30
31 114 analyzed for the concentration of metals using ICP-OES or/and ICP-MS with an external or standard addition
32
33 115 calibration depending on the laboratory procedure. Dilution with 2 % (v/v) HNO_3 was required if the concentration
34
35 116 exceeded the calibration range. The reagent blank was used for background correction.
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60




Figure 2: Tessier sequential extraction procedure flowchart

2.3 Data manipulation

A robust estimator based on the DerSimonian-Laird model was used for the assignment of the consensus value between the various laboratory data sets. NIST Consensus Builder software (<https://consensus.nist.gov/app/nicob>) was used for combining results from two laboratories using mean, standard deviation, and degrees of freedom from each laboratory (30). The DerSimonian-Laird model is often used for establishing an assigned value (*i.e.* key comparison reference value) and associated standard uncertainty for Consultative Committee for Amount of Substance (CCQM) key comparison data (31). The uncertainty values reported here are standard uncertainty at coverage factor $k = 1$ and expanded to cover 95 % confidence interval using k equals to 2. The approach of using DerSimonian-Laird model was also used for the assigned values of extracted elements determined recently for HISS-1, MESS-4 and PACS-3 CRMs by the BCR sequential extraction (21).

The percent distribution of each extraction step and residue was calculated according to Equation 1. The percent distribution fraction is used to evaluate the mobility of elements from the Tessier extraction.

$$\% \text{ Distribution fraction} = \frac{100 \cdot C_{(i,x)}}{\sum (C_{1,x} + C_{2,x} + C_{3,x} + C_{4,x} + C_{\text{residue},x})}$$

Equation 1

In the above equation, $C_{(i,x)}$ is the concentration of analyte (x) from Tessier step i . Denominator in Eq.1 (total combined concentration values) corresponds to the summation of concentrations (C_i) from steps 1 to 4 and residue (C_{residue}).

The standard uncertainty (u_{std} ; $k = 1$) from DerSimonian-Laird model was compared to the target uncertainty (u_{target}) from the Horwitz equation (see Equation 2 below). Horwitz equation is widely used as a target uncertainty in inter-laboratory comparison (32,33), and as criteria for quality control in food and drug analysis (34). The calculation is based on the mass fraction of analyte (C) as follows :

$$\%u_{\text{target}} = 2^{(1 - 0.5 \log C)}$$

Equation 2

The accuracy of Tessier sequential extraction was calculated as a percent recovery (Equation 3) using the combination of the concentrations from four extraction steps and residue, then comparing against the total concentration from the

certified value. If no certified values were available, an average total concentration from the two laboratories was used as reference point.

$$Recovery\ (\%) = \frac{\sum(C_{1,x} + C_{2,x} + C_{3,x} + C_{4,x} + C_{residue,x}) \times 100}{C_{Certified\ value}}$$

Equation 3

3 Results

3.1 Total concentrations of elements in HISS-1, MESS-4, and PACS-3

The concentrations of 17 elements from the three CRMs were analyzed separately by two laboratories and previously reported by Kumkrong et.al (21) as indicated in Table S4 of Electronic Supplementary Material. The percent recovery of each element was calculated against the certified value and presented as a bar graph in Fig. S1 of Electronic Supplementary Material. The total concentrations of elements reported by the two laboratories shows no significant difference between the two digestion methods at 95% confidence interval level by paired T-test (21). The concentration of Cr was significantly low from both digestion methods when a mild digestion procedure was applied (27).

This study reported the concentrations associated with uncertainty (k=2) of some elements that are not certified in the three CRMs including Sb $0.119 \pm 0.049\text{ mg kg}^{-1}$, Sn $0.142 \pm 0.020\text{ mg kg}^{-1}$, and U $0.304 \pm 0.116\text{ mg kg}^{-1}$ for HISS-1, Se $1.18 \pm 0.10\text{ mg kg}^{-1}$ for MESS-4, and Se $1.03 \pm 0.26\text{ mg kg}^{-1}$, and U $2.58 \pm 0.66\text{ mg kg}^{-1}$ for PACS-3.

3.2 Concentrations of elements by Tessier sequential extraction

The concentrations of elements in HISS-1, MESS-4 and PACS-3 from the extraction steps and residue along with the associated uncertainty (k=2) are reported in Table 2 to Table 4. The distribution profiles of the three marine sediment CRMs was calculated for the purpose of interpreting the mobility of the elements and are shown in Fig. 3. The percent distribution was calculated based on the mass fraction of each extraction step and the combined concentration of the four extraction steps and the residue, as mentioned in Equation 1. The concentrations of Cd and Se in HISS-1 were removed from the results shown in Fig. 3 due to an abnormal significantly higher recovery measured for these elements

1
2
3 165 in HISS-1. We note that these two elements in HISS-1 were presented in extremely low concentrations (0.024 mg kg^{-1}
4
5 166 for Cd and 0.050 mg kg^{-1} for Se) (27).
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

167

Table 2 : Assigned value (mg kg⁻¹ dry weight) and associated uncertainty (U, k = 2) of elements in HISS-1

Elements	N(n)*	Tessier # 1		Tessier # 2		Tessier # 3		Tessier # 4		Residue		Combined value **	
		Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U
Li	2(8)	0.105	0.028	0.123	0.035	0.322	0.130	≤ 0.01	-	3.05	2.64	3.60	2.64
Be	2(8)	≤ 0.004	-	≤ 0.004	-	0.004	0.002	0.009	0.002	0.125	0.059	0.125	0.059
V	2(8)	0.133	0.050	≤ 0.01	-	1.23	0.17	0.374	0.126	6.33	2.11	8.06	2.12
Cr	3(13)	≤ 0.003	-	< 0.003	-	1.18	0.21	0.596	0.110	16.8	0.6	18.5	0.6
Mn	2(8)	≤ 0.1	-	1.64	0.44	4.99	0.42	0.502	0.005	59.8	7.1	67.0	7.1
Fe	2(8)	≤ 2	-	6.30	0.15	333	27	18.7	1.0	2,003	110	2,361	113
Ni	2(8)	≤ 0.01	-	0.575	0.107	0.752	0.112	0.345	0.048	1.24	0.29	2.91	0.33
Cu	3(13)	0.113	0.002	0.279	0.069	0.751	0.059	0.711	0.019	0.627	0.311	2.48	0.32
Zn	2(8)	≤ 0.01	-	≤ 0.01	-	1.52	0.14	0.412	0.091	3.14	0.71	5.07	0.72
As	2(8)	0.050	0.026	0.124	0.0003	0.198	0.044	0.027	0.012	0.501	0.403	0.902	0.407
Sr	3(13)	6.97	0.54	52.5	18.2	2.14	0.10	0.266	0.016	21.1	0.5	82.9	18.2
Sn	2(8)	≤ 0.01	-	≤ 0.01	-	≤ 0.01	-	0.054	0.003	0.138	0.114	0.192	0.114
Sb	2(8)	≤ 0.01	-	≤ 0.01	-	0.179	0.023	≤ 0.01	-	0.116	0.005	0.116	0.005

Pb	2(8)	≤ 0.01	-	0.682	0.135	0.622	0.029	0.186	0.069	1.99	0.24	3.48	0.29
U	2(8)	≤ 0.01	-	≤ 0.01	-	0.059	0.022	0.013	0.006	0.157	0.036	0.228	0.043

* N is number of Tessier extraction experiments

n is number of repeated measurements of element concentration

** A combined value is summation of concentrations of an element from T # 1, # 2, # 3, # 4 and residue, and a combined uncertainty is calculated from $u_{std} = 2 \cdot \sqrt{u_{step\ 1}^2 + u_{step\ 2}^2 + u_{step\ 3}^2 + u_{Residue}^2}$ (35)

172

Table 3 : Assigned value (mg kg⁻¹ dry weight) and associated uncertainty (U, k = 2) of elements in MESS-4

Elements	N(n)*	Tessier # 1		Tessier # 2		Tessier # 3		Tessier # 4		Residue		Combined value **	
		Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U
Li	3(11)	0.403	0.086	0.486	0.022	6.80	1.51	2.35	0.28	48.6	7.7	58.7	7.9
Be	2(8)	≤ 0.004	-	0.041	0.024	0.472	0.180	0.070	0.010	1.72	0.13	2.30	0.22
V	2(8)	0.276	0.044	0.471	0.121	21.9	2.3	0.764	0.205	199	25	222	25
Cr	3(11)	0.005	0.001	0.159	0.030	3.55	0.69	1.84	0.40	89.0	10.4	94.5	10.5
Mn	3(11)	2.95	0.26	75.7	7.1	71.5	10.4	11.4	0.24	111	16	273	20
Fe	2(8)	4.26	1.23	327	31	6,152	235	508	108	27,826	6,756	34,818	6,761
Ni	3(11)	0.055	0.009	0.829	0.067	12.3	3.3	5.82	0.89	26.8	4.6	45.8	5.8
Cu	3(11)	0.463	0.015	1.44	0.29	4.42	0.27	7.13	0.93	20.9	2.7	34.4	2.9
Zn	3(11)	0.101	0.044	11.7	0.5	49.4	10.4	13.6	2.1	81.8	8.1	157	13
As	3(11)	0.102	0.014	0.257	0.050	2.57	0.27	0.264	0.060	18.0	3.7	21.2	3.7
Se	2(8)	0.025	0.005	0.205	0.138	0.099	0.090	0.526	0.125	0.219	0.081	1.07	0.27
Sr	3(11)	14.4	0.2	10.7	1.7	10.7	2.0	2.38	0.38	85.4	13.3	124	14
Cd	2(8)	≤ 0.02	-	0.080	0.025	0.049	0.035	0.040	0.072	0.098	0.059	0.241	0.103
Sn	2(8)	≤ 0.01	-	≤ 0.01	-	≤ 0.01	-	≤ 0.01	-	3.23	0.48	3.23	0.48

Sb	2(8)	0.020	0.016	0.028	0.002	0.011	0.007	≤ 0.01	-	1.11	0.18	1.17	0.18
Pb	3(11)	0.014	0.006	2.50	0.34	8.33	1.88	0.271	0.025	10.2	2.7	21.3	3.3
U	2(8)	0.055	0.003	0.237	0.0312	0.353	0.224	0.162	0.004	2.28	0.93	3.09	0.95

* N is number of Tessier extraction experiments

n is number of repeated measurements of element concentration

** A combined value is summation of concentrations of an element from T # 1, # 2, # 3, # 4 and residue, and a combined uncertainty is calculated from $u_{std} = 2 \cdot \sqrt{u_{step\ 1}^2 + u_{step\ 2}^2 + u_{step\ 3}^2 + u_{Residue}^2}$ (35)

176

Table 4 : Assigned value (mg kg⁻¹ dry weight) and associated uncertainty (U, k = 2) of elements in PACS-3

Elements	N(n)*	Tessier # 1		Tessier # 2		Tessier # 3		Tessier # 4		Residue		Combined value **	
		Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U
Li	3(11)	0.518	0.012	0.277	0.083	7.76	0.88	2.37	0.03	19.3	4.9	30.2	5.0
Be	3(11)	≤ 0.004	-	0.032	0.004	0.176	0.050	0.030	0.002	0.750	0.035	0.988	0.061
V	3(11)	0.273	0.310	0.591	0.207	38.3	5.5	4.14	0.48	95.9	19.1	139	20
Cr	3(11)	0.133	0.056	0.917	0.213	20.8	2.3	11.6	1.9	66.1	8.5	99.6	9.0
Mn	3(11)	2.33	0.27	6.55	1.31	49.3	8.0	38.2	1.3	299	39	396	39
Fe	2(8)	4.05	3.31	340	214	7,353	421	4,907	539	24,482	327	37,086	870
Ni	3(11)	0.247	0.080	1.40	0.07	10.6	2.5	4.55	0.73	24.7	4.2	41.4	5.0
Cu	3(11)	4.15	0.83	27.8	7.5	20.6	4.2	223	21	42.2	5.3	318	23
Zn	3(11)	17.2	1.2	140	19	168	38	33.4	3.4	72.0	10.2	430	44
As	3(11)	0.093	0.034	0.363	0.016	10.0	0.2	1.87	0.25	18.5	3.9	30.8	3.9
Se	2(8)	0.024	0.008	0.039	0.011	< 0.02	-	0.888	0.254	0.232	0.005	1.18	0.25
Sr	3(11)	28.1	2.5	11.3	1.9	15.1	1.2	4.08	0.28	174	103	232	103
Cd	2(8)	≤ 0.02	-	1.18	0.33	0.910	0.356	0.085	0.078	0.098	0.060	2.27	0.50
Sn	2(8)	≤ 0.01	-	≤ 0.01	-	1.89	0.10	≤ 0.01	-	24.4	7.9	26.3	7.9

Sb	3(11)	0.046	0.016	0.223	0.009	0.442	0.399	0.079	0.049	13.0	1.3	13.8	1.4
Pb	3(11)	0.919	0.330	37.5	12.3	118	26	22.9	2.4	36.1	8.2	216	30
U	2(8)	≤ 0.01	-	0.584	0.047	0.555	0.211	0.328	0.284	1.17	0.27	2.64	0.45

* N is number of Tessier extraction experiments

n is number of repeated measurements of element concentration

** A combined value is summation of concentrations of an element from T # 1, # 2, # 3, # 4 and residue, and a combined uncertainty is calculated from $u_{std} = 2 \cdot \sqrt{u_{step\ 1}^2 + u_{step\ 2}^2 + u_{step\ 3}^2 + u_{Residue}^2}$ (35)




Figure 3 : Percent distribution (y-scale) of elements (x-scale) by Tessier extraction from A HISS-1, B MESS-4 and C PACS-3

4 Discussion

4.1 Total concentrations of elements in HISS-1, MESS-4, and PACS-3

Fig. 4 plots the total concentrations of elements in the three marine sediment CRMs in comparison to the concentrations in Pelagic clay of the ocean (3) and the earth crust (36). Generally, in decreasing order of abundance, the concentrations of elements were: PACS-3 > MESS-4 > Pelagic clay > Earth crust >> HISS-1. The concentrations of Fe in MESS-4 and PACS-3 were similar and 80 times higher than in HISS-1. HISS-1 contains high amount of silica as quartz and has low accumulation of foreign matters (21,27). PACS-3 showed high concentrations of Zn, Cu, Pb, Sn, As, Sb, Se and Cd, possibly attributable to anthropogenic inputs. MESS-4 showed twice higher concentrations for V and Li when compared to PACS-3.




Figure 4 : Element concentration and elements in HISS-1, MESS-4 and PACS-3 marine sediments, in comparison to Pelagic clay (3) and Earth crust (36)

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

4.2 Distribution profile of elements from Tessier sequential extraction

4.2.1 HISS-1

Considering the released concentrations of elements in HISS-1 (**Table 2**), elements were detected in significant amounts from T # 1 and T # 2, except for Sr. Strontium likely associated with carbonate minerals and dissolved in the acidic solution of T # 2. The concentration of Sr in T # 2 appeared dominant of 63 % and left only 25 % in the residue (**Fig. 3A**). Very little Fe and Mn were extracted in T # 1. High concentrations of Fe and Mn were released in the reducing condition of T # 3; however, majority still remained in the residue (**Fig. 3A**). Additionally, many of the elements were released in T # 3 including Mn (4.95 mg kg⁻¹) > Sr > Zn ≈ V ≈ Cr > Ni ≈ Cu ≈ Pb (0.622 mg kg⁻¹). It is likely that most elements were released with Fe dissolution in reducing condition. Greater than 70 % fraction of Sb, Cr, Mn, Be, Fe, Li, V, and Sn remained in the residues.

4.2.2 MESS-4

The total concentrations of released elements in MESS-4 are 3 to 45 times higher than in HISS-1. The higher concentrations observed in T # 1 were for Sr (14.4 mg kg⁻¹) >> Fe > Mn (2.95 mg kg⁻¹), and trace amounts of Cu (0.463 mg kg⁻¹) > Li > V (0.276 mg kg⁻¹) were observed. At T # 2, significantly high concentrations of Fe and Mn were released, in contrast with HISS-1 where virtually no Fe and Mn were released in T # 1, as noted above. A high fraction of 28 % Mn released at T # 2 suggests that Mn in carbonate minerals was dominant in MESS-4 marine sediment. Besides, high concentrations of Fe (327 mg kg⁻¹) >>> Zn >> Pb > Cu (1.4 mg kg⁻¹) were also observed with carbonate mineral dissolution at T # 2. Other elements that were observed in trace amounts (from 0.1 mg kg⁻¹ to 1 mg kg⁻¹) in T # 2, included As, Cr, Li, Ni, V and U. Several elements were released in T # 3, similarly to HISS-1, but 3 to 116 times larger concentrations were detected in MESS-4 as follows : Fe (6,152 mg kg⁻¹) >>> Mn > Zn > V > Ni ≈ Sr > Pb > Li > Cu > Cr > As (2.57 mg kg⁻¹). At T # 4, apart from high concentrations of Fe, Zn, and Mn, a large fraction (21 %) of the Cu (**Fig. 3B**) was released (7.13 mg kg⁻¹) followed by Ni (5.82 mg kg⁻¹) >> Sr ≈ Li > Cr (1.84 mg kg⁻¹). A significant portion of Sn (100 %), Sb, Cr, V, Li, Fe, Be, and U (74 %) remained in the residue, demonstrating minimal mobility of these elements in MESS-4.

4.2.3 PACS-3

PACS-3 is a highly-contaminated marine sediment with 3 to 200 times higher concentration of the elements analyzed than in HISS-1 (**Fig 4**). At T # 1, elevated concentrations were observed for most elements (except for Be, Cd, Sb, Se, Sn, and U), with particularly high concentrations measured for Zn, Cu, and Pb. At T # 2, higher concentrations of Zn, Pb and Cu were released with carbonate dissolution in comparison to T #1. Similar to HISS-1 and MESS-4, high concentration of Fe (7,353 mg kg⁻¹) was observed at T # 3 together with many elements including Zn (168 mg kg⁻¹) > Pb >> Mn > V > Cr, Cu > Sr > Ni ≈ As (10 mg kg⁻¹). It seems the highest fractions of Pb, Zn, and As are observed at T # 3. Interestingly, Sn was firstly observed in PACS-3 at T # 3 (**Fig. 3C**), suggesting that high accumulation of Sn from the harbor input existed as mono-, di- and tri butyltin in PACS-3 (29); in particular, butyltin is widely used as an antifouling agent in marine paint (37). At T # 4, surprisingly, a high concentration of Fe (4,907 mg kg⁻¹) was observed, indicating Fe associated with sulfide minerals in PACS-3, which is different from HISS-1 and MESS-4. The largest fractions of Cu (70 %) and Se (59 %) were observed at T # 4 (**Fig 3C**), similarly to HISS-1 (**Fig 3A**) and MESS-4 (**Fig 3B**). The most stable elements which remained in the residue are Sb (94 %), Sn, Be, Mn, and Sr (75 %).

4.3 Analytical performance of Tessier sequential extraction on three marine sediment CRMs

4.3.1 Assessment of precision

Since there are no CRM reference data sets readily available on Tessier sequential extraction procedure, assessment of precision and accuracy is needed to benchmark the procedure. The reproducibility from Horwitz formulation (33) was chosen as a target to assess observed precisions in this study. **Fig. 5** illustrates that the majority of standard uncertainty ($u_{\text{std}, k=1}$) of element concentrations determined with the Tessier sequential extractions on HISS-1, MESS-4 and PACS-3 were below the target uncertainty specified by Horwitz equation (Equation 2). Higher than target uncertainties were occasionally observed for lower concentration elements of the three CRMs, which might reflect either stability of the extracts and/or the laboratories ability to perform the measurements at these concentration regimes. Larger uncertainties observed for elements extracted with high concentrations from PACS-3 are unexpected and non-rationalized.

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

257 The distribution of observed relative uncertainties from the four Tessier extraction steps and residue of the three CRMs
258 ranged from 1% to 91% (Cd). The distribution histogram in comparison with target uncertainty from Horwitz model
259 was presented in Fig S2 of Electronic Supplementary Material. Majority (72 %) of uncertainty values in the data set
260 was below the average uncertainty of 13 %, which compares favorably with the average target uncertainty of 16 %
261 (from Horwitz equation, Equation 2).

262
263

For Peer Review




Figure 5 : Standard uncertainty (% , at k=1) of A HISS-1, B MESS-4, and C PACS-3 from Step T # 1, T # 2, T # 3, T # 4 and residue in comparison with target uncertainty from Horwitz equation (red dot).

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

4.3.2 Assessment of accuracy

The assessment of accuracy of the Tessier extraction procedure is non trivial as no quality control materials/ CRMs are available. The accuracy was assessed as recovery by combining the elemental concentrations measured from the four extraction steps and the residue and then comparing against the total concentrations of analytes from either certified values or measured total concentrations obtained here (see Tables S5 to Table S7 of Electronic Supplementary Material). The results of the recovery study are presented as a bar graph in **Fig. 6**. The recovery of HISS-1 was 92 % $\pm$ 40 %, MESS-4 of 101 % $\pm$ 12 % and PACS-3 of 102 % $\pm$ 10 %. A large deviation of recovery values ($\pm$ 40 %) for HISS-1 is due to a low recovery of Cr (62 %), an element which is known to be problematic for analysis of this particular CRM (21,38).




Figure 6 : Percent recovery of combined concentrations of four extraction steps and residue vs elements of A HISS-1, B MESS-4 and C PACS-3 with observed uncertainty (k=2)

5 Conclusions

The four-step Tessier sequential extraction procedure was applied to three widely used marine sediment CRMs and 17 elements were analysed. Different dissolution profiles were observed for HISS-1, MESS-4 and PACS-3. Elements that showed mobility in all the extraction steps included Ni, Cu, Zn, As, Se, Sr, Cd, Pb, and U.

In T # 1, a high amount of Sr was released from the three CRMs, particularly in the highly-contaminated marine sediment PACS-3, which released significant amount of Sr, as well as Zn. In the acidic conditions of T # 2, Sr preferentially released from HISS-1, Mn from MESS-4, and Zn and Cd from PACS-3. Most metals mimic the dissolution pattern of Fe with a significant release under the reducing conditions of extraction step T # 3, particularly Pb, Ni and Zn. Interestingly, in oxidizing conditions (T # 4), Cu and Se are released likely because of their association with iron sulfide, which is being selectively dissolved in this step. In all three sediments, Be, Sb, and Sn remained in the residue. Other elements found in high amounts within the residues of HISS-1 and MESS-4 included Cr, V, Li and Fe.

The precision of Tessier sequential extraction in this study was benchmarked against Horwitz model and was found largely satisfactory. However some elements, which were present in low concentrations, consistently showed high uncertainty. Overall, the percent recovery by the combined element concentrations from the four Tessier extraction steps and residue compared against the certified values shows satisfactory performance of $101 \% \pm 12 \%$ for MESS-4 and $102 \% \pm 10 \%$ for PACS-3, but poor recovery with high uncertainty for HISS-1 of $92 \% \pm 40 \%$. The reference data set produced here using the Tessier sequential extraction procedure carried out on HISS-4, MESS-4 and PACS-3 CRMs will serve as a quality control and method development tool for laboratories.

Acknowledgement

Authors would like to thank the Inorganic Metrology Group and the Environmental Advances in Mining Program at the National Research Council of Canada for support of materials and facilities required to complete this project.

326 Compliance with ethical standards

327 Conflict of interest

328 The authors declare that they have no known competing financial interests or personal relationships that could have
329 appeared to influence the work reported in this paper.

330 References

- 331 1. Van Gestel CAM, Borgman E, Verweij RA, Diez Ortiz M. The influence of soil properties on the toxicity of
332 molybdenum to three species of soil invertebrates. *Ecotoxicol Environ Saf*. 2011;74(1):1–9.
333 <https://doi.org/10.1016/j.ecoenv.2010.10.001>
- 334 2. Moyé J, Picard-Lesteven T, Zouhri L, El Amari K, Hibti M, Benkaddour A. Groundwater assessment and
335 environmental impact in the abandoned mine of Kettara (Morocco). *Environ Pollut*. 2017 ;231:899–907.
336 <https://doi.org/10.1016/j.envpol.2017.07.044>
- 337 3. Li YH, Schoonmaker JE. 9.1 - Chemical Composition and Mineralogy of Marine Sediments. In: Mackenzie
338 FT, editor. *Treatise on Geochemistry: Second Edition*. Elsevier-Pergamon; 2003. p. 1–32.
339 <https://doi.org/10.1016/B978-0-08-095975-7.00701-4>
- 340 4. Worakhunpiset S. Trace Elements in Marine Sediment and Organisms in the Gulf of Thailand. *Int J Environ*
341 *Res Public Health*. 2018;15(4):810. <https://doi.org/10.3390/ijerph15040810>
- 342 5. Joksimović D, Stanković S. Accumulation of trace metals in marine organisms of the southeastern Adriatic
343 coast, Montenegro. *J Serbian Chem Soc*. 2012;77(1):105–17. <https://doi.org/10.2298/JSC110323159J>
- 344 6. Chouvelon T, Strady E, Harmelin-Vivien M, Radakovitch O, Brach-Papa C, Crochet S, et al. Patterns of
345 trace metal bioaccumulation and trophic transfer in a phytoplankton-zooplankton-small pelagic fish marine
346 food web. *Mar Pollut Bull* . 2019;146:1013–30. <https://doi.org/10.1016/j.marpolbul.2019.07.047>
- 347 7. Shefer E, Silverman J, Herut B. Trace metal bioaccumulation in Israeli Mediterranean coastal marine
348 mollusks. *Quat Int*. 2015;390:44–55. <https://doi.org/10.1016/j.quaint.2015.10.030>

- 349 8. Galkus A, Jokšas K, Stakenienė R, Lagunaviciene L. Heavy metal contamination of harbor bottom
350 sediments. *Polish J Environ Stud.* 2012;21(6):1583–94.
- 351 9. Kulkarni P, Chellam S, Ghurye G, Fraser MP. In Situ Generation of Hydrofluoric Acid during Microwave
352 Digestion of Atmospheric Particulate Matter Prior to Trace Element Analysis Using Inductively Coupled
353 Plasma Mass Spectrometry. *Environ Eng Sci.* 2003;20(6):517–31. <https://doi.org/10.1089/109287503770736041>
- 354 10. Lo JM, Sakamoto H. Comparison of the Acid Combinations in Microwave-assisted Digestion of Marine
355 Sediments for Heavy Metal Analyses. *Anal Sci.* 2005 ;21(10):1181–4.
356 <https://doi.org/10.2116/analsci.21.1181>
- 357 11. Salomons W, Förstner U. Trace metal analysis on polluted sediments: Part II: Evaluation of environmental
358 impact. *Environ Technol Lett.* 1980;1(11):506–17. <https://doi.org/10.1080/09593338009384007>
- 359 12. Chester R, Hughes MJ. A chemical technique for the separation of ferro-manganese minerals, carbonate
360 minerals and adsorbed trace elements from pelagic sediments. *Chem Geol.* 1967;2(C):249–62.
361 [https://doi.org/10.1016/0009-2541\(67\)90025-3](https://doi.org/10.1016/0009-2541(67)90025-3)
- 362 13. Barber C. Major and trace element associations in limestones and dolomites. *Chem Geol.* 1974 ;14(4):273–
363 80. [https://doi.org/10.1016/0009-2541\(74\)90064-3](https://doi.org/10.1016/0009-2541(74)90064-3)
- 364 14. Usero J, Gamero M, Morillo J, Gracia I. Comparative study of three sequential extraction procedures for
365 metals in marine sediments. *Environ Int.* 1998;24(4):487–96. [https://doi.org/10.1016/S0160-4120\(98\)00028-](https://doi.org/10.1016/S0160-4120(98)00028-2)
366 2
- 367 15. Gupta SK, Chen KY. Partitioning of trace metals in selective chemical fractions of nearshore sediments.
368 *Environ Lett.* 1975;10(2):129–58. <https://doi.org/10.1080/00139307509435816>
- 369 16. Ure MU, Thomas R, Littlejohn D. Ammonium acetate extracts and their analysis for the speciation of metal
370 ions in soils and sediments. *Int J Environ Anal Chem.* 1993;51(1–4):65–84. <https://doi.org/10.1080/03067319308027612>
- 371
372 17. Quevauviller P, Griepink B, Rauret G. Single and sequential extraction in sediments and soils. *Int J Environ*
373 *Anal Chem.* 1993;51(1–4):231–5. <https://doi.org/10.1080/03067319308027629>

- 1
2
3 374 18. Tessier A, Campbell PGC, Bisson M. Sequential Extraction Procedure for the Speciation of Particulate
4 375 Trace Metals. *Anal Chem*. 1979;51(7):844–51. <https://pubs.acs.org/doi/abs/10.1021/ac50043a017>
5
6
7 376 19. European Commission. The certification of the extractable contents (mass fractions) of Cd, Cr, Cu, Ni, Pb
8 377 and Zn in freshwater sediment following a sequential extraction procedure - BCR 701. Community
9 378 Research. 2001. p. 1–88. [https://op.europa.eu/en/publication-detail/-/publication/02c41803-8f6d-41f2-a73d-](https://op.europa.eu/en/publication-detail/-/publication/02c41803-8f6d-41f2-a73d-c47225ff8c6e)
10 379 [c47225ff8c6e](https://op.europa.eu/en/publication-detail/-/publication/02c41803-8f6d-41f2-a73d-c47225ff8c6e)
11
12
13
14
15 380 20. López-Sánchez JF, Sahuquillo A, Fiedler HD, Rubio R, Rauret G, Muntau H, et al. CRM 601, A stable
16 381 material for its extractable content of heavy metals. *Analyst*. 1998;123(8):1675–7.
17 382 <https://pubs.rsc.org/en/content/articlehtml/1998/an/a802720j>
18
19
20
21 383 21. Kumkrong P, Mercier PHJ, Pihilligawa Gedara I, Mihai O, Tyo DD, Cindy J, et al. Determination of 27
22 384 metals in HISS-1, MESS-4 and PACS-3 marine sediment certified reference materials by the BCR
23 385 sequential extraction. *Talanta*. 2021;221:121543. <https://doi.org/10.1016/j.talanta.2020.121543>
24
25
26
27
28 386 22. Sungur A, Soylak M, Ozcan H. Investigation of heavy metal mobility and availability by the BCR sequential
29 387 extraction procedure: Relationship between soil properties and heavy metals availability. *Chem Speciat*
30 388 *Bioavailab*. 2014;26(4):219–30. <https://doi/abs/10.3184/095422914X14147781158674>
31
32
33
34 389 23. Kubota R, Ohta A, Okai T. Speciation of 38 elements in eight GSJ geochemical sedimentary reference
35 390 materials determined using a sequential extraction technique. *Geochem J*. 2014;48(2):165–88.
36 391 <https://doi.org/10.2343/geochemj.2.0297>
37
38
39
40
41 392 24. Long Y-Y, Hu L-F, Fang C-R, Wu Y-Y, Shen D-S. An evaluation of the modified BCR sequential
42 393 extraction procedure to assess the potential mobility of copper and zinc in MSW. *Microchem J*.
43 394 2009;91(1):1–5. <https://doi.org/10.1016/j.microc.2008.05.006>
44
45
46
47 395 25. Craba L, Brunori C, Galletti M, Cremisini C, Morabito R. Comparison of three sequential extraction
48 396 procedures (original and modified 3 steps BCR procedure) applied to sediments of different origin. *Ann*
49 397 *Chim*. 2004;94(56):409–19. <http://doi.wiley.com/10.1002/adic.200490050>
50
51
52
53 398 26. Tokalioğlu Ş, Kartal Ş, Birol G. Comparison of three sequential extraction procedures for partitioning of
54
55
56
57
58
59
60

- 399 heavy metals in car park dusts. *J Environ Monit.* 2003;5(3):468–76.
- 400 <https://pubs.rsc.org/en/content/articlehtml/2003/em/b300047h>
- 401 27. Berman S, Boyko V, Clancy V, Lam J, Maxwell P, McLaren J, et al. HISS-1 : Marine Sediment Reference
- 402 Material for Trace Metals and other Constituents. 1997. [https://nrc-digital-](https://nrc-digital-repository.canada.ca/eng/view/object/?id=0f2ce620-b76b-4578-ba06-6aae6e632038)
- 403 [repository.canada.ca/eng/view/object/?id=0f2ce620-b76b-4578-ba06-6aae6e632038](https://nrc-digital-repository.canada.ca/eng/view/object/?id=0f2ce620-b76b-4578-ba06-6aae6e632038)
- 404 28. Willie S, Nadeau K, Gedara Pihillagawa I, Yang L, Clancy V, Grinberg P, et al. MESS-4 : Marine Sediment
- 405 Reference Material for Trace Metals and other Constituents. 2014. [https://nrc-digital-](https://nrc-digital-repository.canada.ca/eng/view/object/?id=8a3fd39a-c068-4ce0-820c-d08cf742a20a)
- 406 [repository.canada.ca/eng/view/object/?id=8a3fd39a-c068-4ce0-820c-d08cf742a20a](https://nrc-digital-repository.canada.ca/eng/view/object/?id=8a3fd39a-c068-4ce0-820c-d08cf742a20a)
- 407 29. Willie S, Boyko V, Brophy C, Clancy V, Gedara Pihillagawa I, Grinberg P, et al. PACS-3 : Marine
- 408 Sediment Reference Material for Trace Metals and other Constituents. 2013. [https://nrc-digital-](https://nrc-digital-repository.canada.ca/eng/view/object/?id=6616d284-f101-4b7a-b3c1-a5280e4bd355)
- 409 [repository.canada.ca/eng/view/object/?id=6616d284-f101-4b7a-b3c1-a5280e4bd355](https://nrc-digital-repository.canada.ca/eng/view/object/?id=6616d284-f101-4b7a-b3c1-a5280e4bd355)
- 410 30. Koepke A, Lafarge T, Possolo A, Toman B. NIST Consensus Builder—User’s Manual. NIST, editor. U.S.
- 411 Department of Commerce; 2017. [https://consensus.nist.gov/app_direct/nicob/NISTConsensusBuilder-](https://consensus.nist.gov/app_direct/nicob/NISTConsensusBuilder-UserManual.pdf)
- 412 [UserManual.pdf](https://consensus.nist.gov/app_direct/nicob/NISTConsensusBuilder-UserManual.pdf)
- 413 31. Cox M, Harris P, Ellison S. Use of an “excess-variance” approach for the estimation of a key comparison
- 414 reference value, associated standard uncertainty and degrees of equivalence for CCQM key comparison
- 415 data. *CCQM/11-18.* 2011;1–10.
- 416 <https://pdfs.semanticscholar.org/8f82/448ad8ae426308b1b631748d4e3ed69e891b.pdf>
- 417 32. Beldi G, Jakubowska N, Simoneau C. Report of the interlaboratory comparison organised by the European
- 418 Reference Laboratory for Food Contact Material - Formaldehyde in 3% acetic acid migration solution -
- 419 follow up exercise, EUR 26206 EN,. 2013. [https://op.europa.eu/en/publication-detail/-](https://op.europa.eu/en/publication-detail/-/publication/6e9cdb24-f7e2-41c4-9edf-d21fe9daa9ae/language-en)
- 420 [/publication/6e9cdb24-f7e2-41c4-9edf-d21fe9daa9ae/language-en](https://op.europa.eu/en/publication-detail/-/publication/6e9cdb24-f7e2-41c4-9edf-d21fe9daa9ae/language-en)
- 421 33. Thompson M. The amazing Horwitz function. *amc technical brief.* 2004. Report No.: 17.
- 422 https://www.rsc.org/images/horwitz-function-technical-brief-17_tcm18-214859.pdf
- 423 34. Horwitz W. Evaluation of Analytical Methods Used for Regulation of Foods and Drugs. *Anal Chem.*

- 1
2
3 424 1982;54(1):67–76. <https://pubs.acs.org/doi/abs/10.1021/ac00238a002>
4
5
6 425 35. Ellison SLR, Williams A, editors. Eurachem / CITAC Guide: Quantifying Uncertainty in Analytical
7 426 Measurement. Third edit. Eurachem / CITAC Guide CG 4. Eurachem; 2012. 141 p.
8
9 427 https://www.eurachem.org/images/stories/Guides/pdf/QUAM2012_P1.pdf
10
11
12 428 36. Fleischer M. Recent estimates of the abundance of the elements in the earth's crust. Geol Survey Circ 285.
13
14 429 1953;1–14. <https://pubs.er.usgs.gov/publication/cir285>
15
16
17 430 37. Filipkowska A, Kowalewska G, Pavoni B. Organotin compounds in surface sediments of the Southern
18 431 Baltic coastal zone: A study on the main factors for their accumulation and degradation. Environ Sci Pollut
19 432 Res. 2014;21(3):2077–87. <https://doi.org/10.1007/s11356-013-2115-x>
20
21
22
23 433 38. Kelly WR, Murphy KE, Becker DA, Mann JL. Determination of Cr in certified reference material HISS-1,
24 434 marine sediment, by cold plasma isotope dilution ICP-MS and INAA: Comparison of microwave versus
25 435 closed (Carius) tube digestion. J Anal At Spectrom. 2003;18(2):166–9.
26
27 436 <https://pubs.rsc.org/en/content/articlehtml/2003/ja/b210056h>
28
29
30
31 437

Tessier sequential extraction on 17 elements from three marine sediment certified reference materials (HISS-1, MESS-4 and PACS-3)

Paramee Kumkrong*, Ovi Mihai, Patrick H.J. Mercier, Indu Gedara Pihilligawa, Daniel D. Tyo, and Zoltan Mester

National Research Council Canada, 1200 Montreal Road Ottawa, Ontario Canada, K1A 0R6

*Corresponding author : paramee.kumkrong@nrc-crnc.gc.ca

ORCID : 0000-0002-9281-9974

Abstract

The four-step Tessier sequential extraction procedure is a well-known approach used for environmental and geochemical studies in soil and sediments. However, a lack of reference materials limits its use making implementation and quality control cumbersome. This study applied Tessier sequential extraction to three globally used marine sediment certified reference materials (CRMs) including HISS-1, MESS-4 and PACS-3 with varying levels of contamination.

The study analyzed the distribution of 17 elements throughout the extraction phases. Overall, the percent recovery (sum of steps vs total metal concentration) of all analyzed elements in Tessier extraction was $92 \% \pm 40 \%$ in HISS-1, $101 \% \pm 12 \%$ in MESS-4 and $102 \% \pm 10 \%$ in PACS-3. The observed uncertainty of the individual elemental concentrations averaged of 13 %, which compares favorably with the 16 % target uncertainty derived from Horwitz equation. The reference data set produced here using the Tessier sequential extraction procedure will serve as a quality control and method development tool for laboratories.

Keywords: Distribution. Certified reference material. Sequential extraction. Recovery. Uncertainty.

24	Table of contents	
25	Abstract	1
26	Abbreviation	3
27	1 Introduction	5
28	2. Materials and methods	6
29	2.1 Method for sediment and residue digestion	7
30	2.2 Tessier sequential extraction procedure	8
31	2.3 Data manipulation	10
32	3 Results	11
33	3.1 Total concentrations of elements in HISS-1, MESS-4, and PACS-3	11
34	3.2 Concentrations of elements by Tessier sequential extraction	11
35	4 Discussion	20
36	4.1 Total concentrations of elements in HISS-1, MESS-4, and PACS-3	20
37	4.2 Distribution profile of elements from Tessier sequential extraction	21
38	4.2.1 HISS-1	21
39	4.2.2 MESS-4	21
40	4.2.3 PACS-3	22
41	4.3 Analytical performance of Tessier sequential extraction on three marine sediment CRMs	22
42	4.3.1 Assessment of precision	22
43	4.3.2 Assessment of accuracy	25
44	5 Conclusions	27
45	Acknowledgement	27
46	Compliance with Ethical Standards	28
47	Conflict of interest	28
48	References	28

1
2
3 51 **Abbreviation**
4
5

6		
7		
8	°C	Degree Celsius
9		
10	µm	Micrometer
11		
12	BCR	Community Bureau of Reference
13		
14	CI	Confidence interval
15		
16	CRC	Collision reaction cell
17		
18	CRM(s)	Certified reference material(s)
19		
20	DIW	Deionized water
21		
22	eV	Electron volt
23		
24	Ext sol#	Extraction solution
25		
26	g	Gram
27		
28	h	Hour
29		
30	ICP	Inductively coupled plasma
31		
32	ISO	International Organization for Standardization
33		
34	L	Litre
35		
36	log	Logarithm
37		
38	mol L ⁻¹	Mole per litre
39		
40	mg	Milligram
41		
42	min	Minute
43		
44	mL	Millilitre
45		
46	MS	Mass spectrometry
47		
48	MΩ-cm	Megohm-centimeter
49		
50	na	Not applicable
51		
52	nd	Not detected
53		
54	nr	Not reported
55		
56	NRC	National Research Council Canada
57		

OES	Optical emission spectrometry
rpm	Revolutions per minute
RSD	Relative standard deviation
SD	Standard deviation
S_R	Reproducibility standard deviation
u_{std}	Standard uncertainty
U	Expanded uncertainty
v	Volume
w	Weight

52

53

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

1 Introduction

The association between metals bound to sediments is of interest because their potential release could pose significant environmental challenges (1,2). Marine sediments are mainly composed of biogenous, hydrogenous, and lithogenous sediments (3), and also receive foreign materials from inland activities such as weathering and anthropogenic inputs. A high accumulation of toxic metals in marine sediments is widely observed in areas of high human activity, for example harbors (4–8).

In general, assessment of sediment quality is achieved by quantitative analysis of total toxic metal content (9–11). However, not all measured total metal content in sediments is actually posing an environmental hazard, so better approaches that estimate the labile metal fraction and the conditions of their release are needed. The importance of such investigations brought the development of methodologies based on fractionating toxic metals from solid phases by partial dissolution / desorption / extraction (12–17).

Among sequential extraction procedures, the Tessier and the Community Bureau of Reference (BCR) sequential extractions are most commonly employed (18,19). The publications citing Tessier and BCR sequential extraction in Scopus database (<https://www.scopus.com>, accessed on June 7, 2020) showed results of 3,625 and 2,332, respectively. Tessier expressed significant growth until it reached a plateau in 2011, however, BCR extraction method has continued gaining popularity (**Fig.1**).

Tessier sequential extraction procedure was published in 1979 for the determination of Cd, Cr, Cu, Fe, Mn, Ni, Pb, and Zn in river sediments from Quebec, Canada (18). The BCR method was developed with the objective to harmonize the sequential extraction procedure and be used in conjunction with the production of BCR 601 and BCR 701 CRMs (river sediment reference materials) that certified Cd, Cr, Cu, Ni, Pb, and Zn (19,20). There is a plethora of data available from the BCR on method characteristics (21–26) but the availability of such information is much limited for the Tessier extraction approach using reference materials.




Figure 1 : The number of publications of Tessier and BCR sequential extraction procedure since 1979 to 2020 (Source: Scopus database, searched under keywords: “sequential extraction BCR” and “sequential extraction Tessier” (<https://www.scopus.com>, accessed on June 07, 2020))

Considering the lack of CRMs and reference data on precision and accuracy upon application of a given extraction method, we have performed the Tessier extraction procedure on 17 elements (As, Be, Cd, Cr, Cu, Fe, Li, Mn, Ni, Pb, Sb, Se, Sn, Sr, U, V, and Zn) in three marine sediment CRMs produced by National Research Council Canada (NRC). The three CRMs included: HISS-1, a clean sediment; MESS-4, a moderately-contaminated sediment; and PACS-3, a highly-contaminated marine sediment (27–29). These were certified for total concentrations of numerous elements and have been used as quality control and for method validation purpose in analytical laboratories worldwide.

2. Materials and methods

HISS-1 CRM was sampled from Hibernia shelf, off the coast of Newfoundland, Canada (27), MESS-4 from Beaufort Sea, Canadian Arctic (28) and PACS-3 from Esquimalt harbor, British Columbia, Canada (29). The CRMs were sieved to particle sizes of less than 125 μm . Homogeneity and stability was assessed. HISS-1 CRM has served as marine sediment reference materials since 1997, PACS-3 since 2013 and MESS-4 since 2014. The three marine sediment CRMs were used “as is” and the moisture content was corrected to obtain concentration based on dry weight.

Laboratory ware, chemicals, and instruments are listed in Table S1 of [Supplementary MaterialElectronic Supplementary Material](#). The extraction solution (Ext sol #) and its preparation in each step are provided in **Table 1**.

Table 1: Extraction solutions for Tessier sequential extraction

Ext sol #	Reagent	Reagent preparation
1	1 mol L ⁻¹ sodium acetate solution (NaOAc), pH 8.2	Dissolve 41 g of NaOAc in 450 mL of deionized water (DIW), adjust to pH 8.2 using acetic acid (AcOH) and dilute to 500 mL with DIW
2	1 mol L ⁻¹ NaOAc, pH 5.0	Same as Ext sol# 1 but adjust to pH 5 using AcOH
3	0.04 mol L ⁻¹ hydroxylammonium chloride solution (NH ₂ OH.HCl) in 25 % (v/v) AcOH	Dissolve 1.4 g of NH ₂ OH.HCl in 500 mL of 25 % (v/v) AcOH
4	0.02 mol L ⁻¹ nitric acid (HNO ₃)	- Add 6.3 mL of 69 % (w/w) HNO ₃ and make up to 1,000 mL with DIW to obtain 0.1 mol L ⁻¹ HNO ₃ - Pipette 20 mL of 0.1 mol L ⁻¹ HNO ₃ and make up to 100 mL with DIW to obtain 0.02 mol L ⁻¹ HNO ₃
	30 % (w/v) hydrogen peroxide (H ₂ O ₂) at pH 2	Add few drops of 69% (w/w) HNO ₃ to 30% (w/v) H ₂ O ₂ and check pH with pH strip to obtain pH 2
	3.2 mol L ⁻¹ ammonium acetate solution (NH ₄ OAc) in 20 % (v/v) HNO ₃	Dissolve 24.7 g of NH ₄ OAc in 100 mL of 20 % (v/v) of HNO ₃

2.1 Method for sediment and residue digestion

The following digestion method is used for total sediment digestion of the raw material CRMs and reported results in Section 3.1 Total concentrations of elements and Section 3.2 Concentrations of elements by Tessier sequential extraction in Table 2 to 4 for the residues. Two laboratories ([Lab A](#) and [Lab B](#)) used their routine method for digestion. ~~Lab A and Lab B used the microwaved digestion method~~ as described by Kumkrong

et.al. (21). In brief Lab A used 0.25 g of sample and a mixture of HNO₃, hydrochloric acid (HCl) and hydrogen fluoride (HF) ~~rammed~~. The excess amount of HF was destroyed by adding boric acid (H₃BO₃)₂. ~~After and re-digested in the~~ microwave digestion, the samples were placed in a hot block, ~~evaporated at~~ to obtain 5 mL volume and diluted to 50 g using 2 % (v/v) HNO₃.

~~Lab B modified the U.S. EPA. Method 3052 (31) as Lab B used follows: 0.25 g of sample or 1.0 g of residue was weighed in pre-cleaned Teflon microwave digestion vessels, then 0.25 g of raw material or 1 g of residue and digested using a mixture of 10 mL of 69 % (w/w) HNO₃, 3 mL of 48 % (w/v) HF, and 3 mL of 30 % (w/w) H₂O₂ were added and the vessels were placed in a fume hood overnight. The next day, the vessels were heated on a hot plate at 80 °C for 24 h. 1 mL of 37 % (w/w) HCl was added and digested using the microwave programmed at 5 min ramping to 360 Watts and hold for 5 min, another 5 min ramping to 850 Watts and hold for 20 min, then 0 Watt for the 35 min cool-down cycle. The digested solution was transferred to a Teflon tube. The opened digestion was performed after microwave digestion using a mixture of HNO₃ and HCl on a hot block and evaporated to approximately near 1 mL. Another 2 mL of 69 % (w/w) HNO₃ and 3 mL of 37 % (w/w) HCl were added to the tube and heated at 80 °C for 60 min. Finally, the and the digested samples were~~ diluted to 100 g using DIW.

The elemental analysis of digestrates was conducted using inductively coupled plasma-optical emission spectrometry (ICP-OES) and inductively coupled plasma-mass spectrometry (ICP-MS). The conditions of ICP-OES and ICP-MS are stated in the Electronic Supplementary Material Table S2.

2.2 Tessier sequential extraction procedure

The extraction procedure is presented in **Fig 2.** and detailed in Table S32 of ~~Supplementary Material~~ Electronic Supplementary Material. The extraction tube was weighed before and after adding a new extraction solution for a weight correction. The extraction steps were labeled as Tessier step 1 - exchangeable (T # 1), step 2 - bound to carbonates (T # 2), step 3 - bound to Fe-Mn oxides (T # 3), step 4 - bound to organic matter and sulfide minerals (T # 4), and residue. The fractions were analyzed for the concentration of metals using ~~an inductively coupled plasma optical emission spectrometer (ICP-OES) or/and inductively coupled plasma mass spectrometer (ICP-MS)~~ with an external or standard addition calibration depending on the laboratory procedure. Dilution with 2 % (v/v) HNO₃ was required if the concentration exceeded the calibration range. The reagent blank was used for background correction.




Figure 2: Tessier sequential extraction procedure flowchart

2.3 Data manipulation

A robust estimator based on the DerSimonian-Laird model was used for the assignment of the consensus value between the various laboratory data sets. NIST Consensus Builder software (<https://consensus.nist.gov/app/nicob>) was used for combining results from two laboratories using mean, standard deviation, and degrees of freedom from each laboratory (30). The DerSimonian-Laird model is often used for establishing an assigned value (*i.e.* key comparison reference value) and associated standard uncertainty for Consultative Committee for Amount of Substance (CCQM) key comparison data (31). The uncertainty values reported here are standard uncertainty at coverage factor $k = 1$ and expanded to cover 95 % confidence interval using k equals to 2. The approach of using DerSimonian-Laird model was also used for the assigned values of extracted elements determined recently for HISS-1, MESS-4 and PACS-3 CRMs by the BCR sequential extraction (21).

The percent distribution of each extraction step and residue was calculated according to Equation 1. The percent distribution fraction is used to evaluate the mobility of elements from the Tessier extraction.

$$\% \text{ Distribution fraction} = \frac{100 \cdot C_{(i,x)}}{\sum (C_{1,x} + C_{2,x} + C_{3,x} + C_{4,x} + C_{\text{residue}, x})}$$

Equation 1

In the above equation, $C_{(i,x)}$ is the concentration of analyte (x) from Tessier step i . Denominator in Eq.1 (total combined concentration values) corresponds to the summation of concentrations (C_i) from steps 1 to 4 and residue (C_{residue}).

The standard uncertainty (u_{std} ; $k = 1$) from DerSimonian-Laird model was compared to the target uncertainty (u_{target}) from the Horwitz equation (see Equation 2 below). Horwitz equation is widely used as a target uncertainty in inter-laboratory comparison (32,33), and as criteria for quality control in food and drug analysis (34). The calculation is based on the mass fraction of analyte (C) as follows :

$$\%u_{\text{target}} = 2^{(1 - 0.5 \log C)}$$

Equation 2

The accuracy of Tessier sequential extraction was calculated as a percent recovery (Equation 3) using the combination of the concentrations from four extraction steps and residue, then comparing against the total concentration from the

certified value. If no certified values were available, an average total concentration from the two laboratories was used as reference point.

$$Recovery\ (\%) = \frac{\sum(C_{1,x} + C_{2,x} + C_{3,x} + C_{4,x} + C_{residue,x}) \times 100}{C_{Certified\ value}}$$

Equation 3

3 Results

3.1 Total concentrations of elements in HISS-1, MESS-4, and PACS-3

The concentrations of 17 elements from the three CRMs were analyzed ~~by two digestions performed separately by each of the two laboratories as previously and reported by Kumkrong et.al (21) and indicated are reported in~~ Table S43 of ~~Supplementary Material~~Electronic Supplementary Material. The percent recovery of each element was calculated against the certified value and presented as a bar graph in Fig. S1 of ~~Supplementary Material~~Electronic Supplementary Material. ~~Paired t-test was used to access the equivalence of the results and showed~~The total concentrations of elements reported by the two laboratories shows no significant difference between the two digestion methods; ~~and between each digestion and certified values at 95% confidence interval level by paired T-test (Table S4 of Supplementary Material)(21). The concentration of Cr was significantly low from both digestion methods. The low concentration of Cr in HISS-1 is reported~~ when a mild digestion procedure was applied- (27).

This study reported the concentrations associated with uncertainty (k=2) of some elements that are not certified in the three CRMs including Sb 0.119 ± 0.049 mg kg⁻¹, Sn 0.142 ± 0.020 mg kg⁻¹, and U 0.304 ± 0.116 mg kg⁻¹ for HISS-1, Se 1.18 ± 0.10 mg kg⁻¹ for MESS-4, and Se 1.03 ± 0.26 mg kg⁻¹, and U 2.58 ± 0.66 mg kg⁻¹ for PACS-3.

3.2 Concentrations of elements by Tessier sequential extraction

The concentrations of elements in HISS-1, MESS-4 and PACS-3 from the extraction steps and residue along with the associated uncertainty (k=2) are reported in Table 2 to Table 4. The distribution profiles of the three marine sediment CRMs was calculated for the purpose of interpreting the mobility of the elements and are shown in Fig. 3. The percent distribution was calculated based on the mass fraction of each extraction step and the combined concentration of the

1
2
3 175 four extraction steps and the residue, as mentioned in Equation 1. The concentrations of Cd and Se in HISS-1 were
4
5 176 removed from the results shown in **Fig. 3** due to an abnormal significantly higher recovery measured for these elements
6
7 177 in HISS-1. We note that these two elements in HISS-1 were presented in extremely low concentrations (0.024 mg/
8
9 178 kg^{-1} for Cd and 0.050 mg/ kg^{-1} for Se) (27).
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

179

Table 2 : Assigned value (mg kg⁻¹ dry weight) and associated uncertainty (U, k = 2) of elements in HISS-1

Elements	N(n)*	Tessier # 1		Tessier # 2		Tessier # 3		Tessier # 4		Residue		Combined value **	
		Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U
Li	2(8)	0.105	0.028	0.123	0.035	0.322	0.130	≤ 0.01	-	3.05	2.64	3.60	2.64
Be	2(8)	≤ 0.004	-	≤ 0.004	-	0.004	0.002	0.009	0.002	0.125	0.059	0.125	0.059
V	2(8)	0.133	0.050	≤ 0.01	-	1.23	0.17	0.374	0.126	6.33	2.11	8.06	2.12
Cr	3(13)	≤ 0.003	-	< 0.003	-	1.18	0.21	0.596	0.110	16.8	0.6	18.5	0.6
Mn	2(8)	≤ 0.1	-	1.64	0.44	4.99	0.42	0.502	0.005	59.8	7.1	67.0	7.1
Fe	2(8)	≤ 2	-	6.30	0.15	333	27	18.7	1.0	2,003	110	2,361	113
Ni	2(8)	≤ 0.01	-	0.575	0.107	0.752	0.112	0.345	0.048	1.24	0.29	2.91	0.33
Cu	3(13)	0.113	0.002	0.279	0.069	0.751	0.059	0.711	0.019	0.627	0.311	2.48	0.32
Zn	2(8)	≤ 0.01	-	≤ 0.01	-	1.52	0.14	0.412	0.091	3.14	0.71	5.07	0.72
As	2(8)	0.050	0.026	0.124	0.0003	0.198	0.044	0.027	0.012	0.501	0.403	0.902	0.407
Sr	3(13)	6.97	0.54	52.5	18.2	2.14	0.10	0.266	0.016	21.1	0.5	82.9	18.2
Sn	2(8)	≤ 0.01	-	≤ 0.01	-	≤ 0.01	-	0.054	0.003	0.138	0.114	0.192	0.114
Sb	2(8)	≤ 0.01	-	≤ 0.01	-	0.179	0.023	≤ 0.01	-	0.116	0.005	0.116	0.005

Pb	2(8)	≤ 0.01	-	0.682	0.135	0.622	0.029	0.186	0.069	1.99	0.24	3.48	0.29
U	2(8)	≤ 0.01	-	≤ 0.01	-	0.059	0.022	0.013	0.006	0.157	0.036	0.228	0.043

* N is number of Tessier extraction experiments

n is number of repeated measurements of element concentration

** A combined value is summation of concentrations of an element from T # 1, # 2, # 3, # 4 and residue, and a combined uncertainty is calculated from $u_{std} = 2 \cdot \sqrt{u_{step\ 1}^2 + u_{step\ 2}^2 + u_{step\ 3}^2 + u_{Residue}^2}$ (35)

184

Table 3 : Assigned value (mg kg⁻¹ dry weight) and associated uncertainty (U, k = 2) of elements in MESS-4

Elements	N(n)*	Tessier # 1		Tessier # 2		Tessier # 3		Tessier # 4		Residue		Combined value **	
		Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U
Li	3(11)	0.403	0.086	0.486	0.022	6.80	1.51	2.35	0.28	48.6	7.7	58.7	7.9
Be	2(8)	≤ 0.004	-	0.041	0.024	0.472	0.180	0.070	0.010	1.72	0.13	2.30	0.22
V	2(8)	0.276	0.044	0.471	0.121	21.9	2.3	0.764	0.205	199	25	222	25
Cr	3(11)	0.005	0.001	0.159	0.030	3.55	0.69	1.84	0.40	89.0	10.4	94.5	10.5
Mn	3(11)	2.95	0.26	75.7	7.1	71.5	10.4	11.4	0.24	111	16	273	20
Fe	2(8)	4.26	1.23	327	31	6,152	235	508	108	27,826	6,756	34,818	6,761
Ni	3(11)	0.055	0.009	0.829	0.067	12.3	3.3	5.82	0.89	26.8	4.6	45.8	5.8
Cu	3(11)	0.463	0.015	1.44	0.29	4.42	0.27	7.13	0.93	20.9	2.7	34.4	2.9
Zn	3(11)	0.101	0.044	11.7	0.5	49.4	10.4	13.6	2.1	81.8	8.1	157	13
As	3(11)	0.102	0.014	0.257	0.050	2.57	0.27	0.264	0.060	18.0	3.7	21.2	3.7
Se	2(8)	0.025	0.005	0.205	0.138	0.099	0.090	0.526	0.125	0.219	0.081	1.07	0.27
Sr	3(11)	14.4	0.2	10.7	1.7	10.7	2.0	2.38	0.38	85.4	13.3	124	14
Cd	2(8)	≤ 0.02	-	0.080	0.025	0.049	0.035	0.040	0.072	0.098	0.059	0.241	0.103
Sn	2(8)	≤ 0.01	-	≤ 0.01	-	≤ 0.01	-	≤ 0.01	-	3.23	0.48	3.23	0.48

Sb	2(8)	0.020	0.016	0.028	0.002	0.011	0.007	≤ 0.01	-	1.11	0.18	1.17	0.18
Pb	3(11)	0.014	0.006	2.50	0.34	8.33	1.88	0.271	0.025	10.2	2.7	21.3	3.3
U	2(8)	0.055	0.003	0.237	0.0312	0.353	0.224	0.162	0.004	2.28	0.93	3.09	0.95

* N is number of Tessier extraction experiments

n is number of repeated measurements of element concentration

** A combined value is summation of concentrations of an element from T # 1, # 2, # 3, # 4 and residue, and a combined uncertainty is calculated from $u_{std} = 2 \cdot \sqrt{u_{step\ 1}^2 + u_{step\ 2}^2 + u_{step\ 3}^2 + u_{Residue}^2}$ (35)

188

Table 4 : Assigned value (mg kg⁻¹ dry weight) and associated uncertainty (U, k = 2) of elements in PACS-3

Elements	N(n)*	Tessier # 1		Tessier # 2		Tessier # 3		Tessier # 4		Residue		Combined value **	
		Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U	Assigned value	U
Li	3(11)	0.518	0.012	0.277	0.083	7.76	0.88	2.37	0.03	19.3	4.9	30.2	5.0
Be	3(11)	≤ 0.004	-	0.032	0.004	0.176	0.050	0.030	0.002	0.750	0.035	0.988	0.061
V	3(11)	0.273	0.310	0.591	0.207	38.3	5.5	4.14	0.48	95.9	19.1	139	20
Cr	3(11)	0.133	0.056	0.917	0.213	20.8	2.3	11.6	1.9	66.1	8.5	99.6	9.0
Mn	3(11)	2.33	0.27	6.55	1.31	49.3	8.0	38.2	1.3	299	39	396	39
Fe	2(8)	4.05	3.31	340	214	7,353	421	4,907	539	24,482	327	37,086	870
Ni	3(11)	0.247	0.080	1.40	0.07	10.6	2.5	4.55	0.73	24.7	4.2	41.4	5.0
Cu	3(11)	4.15	0.83	27.8	7.5	20.6	4.2	223	21	42.2	5.3	318	23
Zn	3(11)	17.2	1.2	140	19	168	38	33.4	3.4	72.0	10.2	430	44
As	3(11)	0.093	0.034	0.363	0.016	10.0	0.2	1.87	0.25	18.5	3.9	30.8	3.9
Se	2(8)	0.024	0.008	0.039	0.011	< 0.02	-	0.888	0.254	0.232	0.005	1.18	0.25
Sr	3(11)	28.1	2.5	11.3	1.9	15.1	1.2	4.08	0.28	174	103	232	103
Cd	2(8)	≤ 0.02	-	1.18	0.33	0.910	0.356	0.085	0.078	0.098	0.060	2.27	0.50
Sn	2(8)	≤ 0.01	-	≤ 0.01	-	1.89	0.10	≤ 0.01	-	24.4	7.9	26.3	7.9

Sb	3(11)	0.046	0.016	0.223	0.009	0.442	0.399	0.079	0.049	13.0	1.3	13.8	1.4
Pb	3(11)	0.919	0.330	37.5	12.3	118	26	22.9	2.4	36.1	8.2	216	30
U	2(8)	≤ 0.01	-	0.584	0.047	0.555	0.211	0.328	0.284	1.17	0.27	2.64	0.45

* N is number of Tessier extraction experiments

n is number of repeated measurements of element concentration

** A combined value is summation of concentrations of an element from T # 1, # 2, # 3, # 4 and residue, and a combined uncertainty is calculated from $u_{std} = 2 \cdot \sqrt{u_{step\ 1}^2 + u_{step\ 2}^2 + u_{step\ 3}^2 + u_{Residue}^2}$ (35)




Figure 3 : Percent distribution (y-scale) of elements (x-scale) by Tessier extraction from A HISS-1, B MESS-4 and C PACS-3

4 Discussion

4.1 Total concentrations of elements in HISS-1, MESS-4, and PACS-3

Fig. 4 plots the total concentrations of elements in the three marine sediment CRMs in comparison to the concentrations in Pelagic clay of the ocean (3) and the earth crust (36). Generally, in decreasing order of abundance, the concentrations of elements were: PACS-3 > MESS-4 > Pelagic clay > Earth crust >> HISS-1. The concentrations of Fe in MESS-4 and PACS-3 were similar and 80 times higher than in HISS-1. HISS-1 contains high amount of silica as quartz and has low accumulation of foreign matters (21,27). PACS-3 showed high concentrations of Zn, Cu, Pb, Sn, As, Sb, Se and Cd, possibly attributable to anthropogenic inputs. MESS-4 showed twice higher concentrations for V and Li when compared to PACS-3.




Figure 4 : Element concentration and elements in profile of HISS-1, MESS-4 and PACS-3 marine sediments, in comparison to Pelagic clay (3) and Earth crust (36)

4.2 Distribution profile of elements from Tessier sequential extraction

4.2.1 HISS-1

Considering the released concentrations of elements in HISS-1 (**Table 2**), elements were detected in significant amounts from T # 1 and T # 2, except for Sr. Strontium likely associated with carbonate minerals and dissolved in the acidic solution of T # 2. The concentration of Sr in T # 2 appeared dominant of 63 % and left only 25 % in the residue (**Fig. 3A**). Very little Fe and Mn were extracted in T # 1. High concentrations of Fe and Mn were released in the reducing condition of T # 3; however, majority still remained in the residue (**Fig. 3A**). Additionally, many of the elements were released in T # 3 including Mn (4.95 mg kg^{-1}) > Sr > Zn $\approx$ V $\approx$ Cr > Ni $\approx$ Cu $\approx$ Pb (0.622 mg kg^{-1}). It is likely that most elements were released with Fe dissolution in reducing condition. Greater than 70 % fraction of Sb, Cr, Mn, Be, Fe, Li, V, and Sn remained in the residues.

4.2.2 MESS-4

The total concentrations of released elements in MESS-4 are 3 to 45 times higher than in HISS-1. The higher concentrations observed in T # 1 were for Sr (14.4 mg kg^{-1}) >> Fe > Mn (2.95 mg kg^{-1}), and trace amounts of Cu (0.463 mg kg^{-1}) > Li > V (0.276 mg kg^{-1}) were observed. At T # 2, significantly high concentrations of Fe and Mn were released, in contrast with HISS-1 where virtually no Fe and Mn were released in T # 1, as noted above. A high fraction of 28 % Mn released at T # 2 suggests that Mn in carbonate minerals was dominant in MESS-4 marine sediment. Besides, high concentrations of Fe (327 mg kg^{-1}) >>> Zn >> Pb > Cu (1.4 mg kg^{-1}) were also observed with carbonate mineral dissolution at T # 2. Other elements that were observed in trace amounts (from 0.1 mg kg^{-1} to 1 mg kg^{-1}) in T # 2, included As, Cr, Li, Ni, V and U. Several elements were released in T # 3, similarly to HISS-1, but 3 to 116 times larger concentrations were detected in MESS-4 as follows : Fe ($6,152 \text{ mg kg}^{-1}$) >>> Mn > Zn > V > Ni $\approx$ Sr > Pb > Li > Cu > Cr > As (2.57 mg kg^{-1}). At T # 4, apart from high concentrations of Fe, Zn, and Mn, a large fraction (21 %) of the Cu (**Fig. 3B**) was released (7.13 mg kg^{-1}) followed by Ni (5.82 mg kg^{-1}) >> Sr $\approx$ Li > Cr (1.84 mg kg^{-1}). A significant portion of Sn (100 %), Sb, Cr, V, Li, Fe, Be, and U (74 %) remained in the residue, demonstrating minimal mobility of these elements in MESS-4.

4.2.3 PACS-3

PACS-3 is a highly-contaminated marine sediment with 3 to 200 times higher concentration of the elements analyzed than in HISS-1 (**Fig 4**). At T # 1, elevated concentrations were observed for most elements (except for Be, Cd, Sb, Se, Sn, and U), with particularly high concentrations measured for Zn, Cu, and Pb. At T # 2, higher concentrations of Zn, Pb and Cu were released with carbonate dissolution in comparison to T #1. Similar to HISS-1 and MESS-4, high concentration of Fe (7,353 mg kg⁻¹) was observed at T # 3 together with many elements including Zn (168 mg kg⁻¹) > Pb >> Mn > V > Cr, Cu > Sr > Ni ≈ As (10 mg kg⁻¹). It seems the highest fractions of Pb, Zn, and As are observed at T # 3. Interestingly, Sn was firstly observed in PACS-3 at T # 3 (**Fig. 3C**), suggesting that high accumulation of Sn from the harbor input existed as mono-, di- and tri butyltin in PACS-3 (29); in particular, butyltin is widely used as an antifouling agent in marine paint (37). At T # 4, surprisingly, a high concentration of Fe (4,907 mg kg⁻¹) was observed, indicating Fe associated with sulfide minerals in PACS-3, which is different from HISS-1 and MESS-4. The largest fractions of Cu (70 %) and Se (59 %) were observed at T # 4 (**Fig 3C**), similarly to HISS-1 (**Fig 3A**) and MESS-4 (**Fig 3B**). The most stable elements which remained in the residue are Sb (94 %), Sn, Be, Mn, and Sr (75 %).

4.3 Analytical performance of Tessier sequential extraction on three marine sediment CRMs

4.3.1 Assessment of precision

Since there are no CRM reference data sets readily available on Tessier sequential extraction procedure, assessment of precision and accuracy is needed to benchmark the procedure. The reproducibility from Horwitz formulation (33) was chosen as a target to assess observed precisions in this study. **Fig. 5** illustrates that the majority of standard uncertainty ($u_{\text{std}, k=1}$) of element concentrations determined with the Tessier sequential extractions on HISS-1, MESS-4 and PACS-3 were below the target uncertainty specified by Horwitz equation (Equation 2). Higher than target uncertainties were occasionally observed for lower concentration elements of the three CRMs, which might reflect either stability of the extracts and/or the laboratories ability to perform the measurements at these concentration regimes. Larger uncertainties observed for elements extracted with high concentrations from PACS-3 are unexpected and non-rationalized.

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

269 The distribution of observed relative uncertainties from the four Tessier extraction steps and residue of the three CRMs
270 ranged from 1% to 91% (Cd). The distribution histogram in comparison with target uncertainty from Horwitz model
271 was presented in Fig S2 of ~~Supplementary Material~~[Electronic Supplementary Material](#). Majority (72 %) of uncertainty
272 values in the data set was below the average uncertainty of 13 %, which compares favorably with the average target
273 uncertainty of 16 % (from Horwitz equation, Equation 2).

274

275

For Peer Review




Figure 5 : Standard uncertainty (% , at k=1) of A HISS-1, B MESS-4, and C PACS-3 from Step T # 1, T # 2, T # 3, T # 4 and residue in comparison with target uncertainty from Horwitz equation (red dot).

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

278 4.3.2 Assessment of accuracy

279 The assessment of accuracy of the Tessier extraction procedure is non trivial as no quality control materials/ CRMs
280 are available. The accuracy was assessed as recovery by combining the elemental concentrations measured from the
281 four extraction steps and the residue and then comparing against the total concentrations of analytes from either
282 certified values or measured total concentrations obtained here (see Tables [S5 to Table S7](#) of [Supplementary](#)
283 [MaterialElectronic Supplementary Material](#)). The results of the recovery study are presented as a bar graph in **Fig. 6**.
284 The recovery of HISS-1 was 92 % $\pm$ 40 %, MESS-4 of 101 % $\pm$ 12 % and PACS-3 of 102 % $\pm$ 10 %. A large deviation
285 of recovery values ($\pm$ 40 %) for HISS-1 is due to a low recovery of Cr (62 %), an element which is known to be
286 problematic for analysis of this particular CRM (21,38).

287




Figure 6 : Percent recovery of combined concentrations of four extraction steps and residue vs elements of A HISS-1, B MESS-4 and C PACS-3 with observed uncertainty ($k=2$)

316 5 Conclusions

317 The four-step Tessier sequential extraction procedure was applied to three widely used marine sediment CRMs and
318 17 elements were analysed. Different dissolution profiles were observed for HISS-1, MESS-4 and PACS-3. Elements
319 that showed mobility in all the extraction steps included Ni, Cu, Zn, As, Se, Sr, Cd, Pb, and U.

320 In T # 1, a high amount of Sr was released from the three CRMs, particularly in the highly-contaminated marine
321 sediment PACS-3, which released significant amount of Sr, as well as Zn. In the acidic conditions of T # 2, Sr
322 preferentially released from HISS-1, Mn from MESS-4, and Zn and Cd from PACS-3. Most metals mimic the
323 dissolution pattern of Fe with a significant release under the reducing conditions of extraction step T # 3, particularly
324 Pb, Ni and Zn. Interestingly, in oxidizing conditions (T # 4), Cu and Se are released likely because of their association
325 with iron sulfide, which is being selectively dissolved in this step. In all three sediments, Be, Sb, and Sn remained in
326 the residue. Other elements found in high amounts within the residues of HISS-1 and MESS-4 included Cr, V, Li and
327 Fe.

328 The precision of Tessier sequential extraction in this study was benchmarked against Horwitz model and was found
329 largely satisfactory. However some elements, which were present in low concentrations, consistently showed high
330 uncertainty. Overall, the percent recovery by the combined element concentrations from the four Tessier extraction
331 steps and residue compared against the certified values shows satisfactory performance of $101 \% \pm 12 \%$ for MESS-
332 4 and $102 \% \pm 10 \%$ for PACS-3, but poor recovery with high uncertainty for HISS-1 of $92 \% \pm 40 \%$. The reference
333 data set produced here using the Tessier sequential extraction procedure carried out on HISS-4, MESS-4 and PACS-
334 3 CRMs will serve as a quality control and method development tool for laboratories.

335 Acknowledgement

336 Authors would like to thank the Inorganic Metrology Group and the Environmental Advances in Mining Program at
337 the National Research Council of Canada for support of materials and facilities required to complete this project.

338 Compliance with ethical standards

339 **Conflict of interest**

340 The authors declare that they have no known competing financial interests or personal relationships that could have
341 appeared to influence the work reported in this paper.

342 **References**

- 343 1. Van Gestel CAM, Borgman E, Verweij RA, Diez Ortiz M. The influence of soil properties on the toxicity of
344 molybdenum to three species of soil invertebrates. *Ecotoxicol Environ Saf*. 2011 ~~Jan 1~~;74(1):1–9.
345 <https://doi.org/10.1016/j.ecoenv.2010.10.001>.
- 346 2. Moyé J, Picard-Lesteven T, Zouhri L, El Amari K, Hibti M, Benkaddour A. Groundwater assessment and
347 environmental impact in the abandoned mine of Kettara (Morocco). *Environ Pollut*. 2017 ~~Dec 1~~;231:899–
348 907. <https://doi.org/10.1016/j.envpol.2017.07.044>.
- 349 3. Li YH, Schoonmaker JE. 9.1 - Chemical Composition and Mineralogy of Marine Sediments. In: Mackenzie
350 FT, editor. *Treatise on Geochemistry: Second Edition*. Elsevier-Pergamon; 2003. p. 1–32.
351 <https://doi.org/10.1016/B978-0-08-095975-7.00701-4>.
- 352 4. Worakhunpiset S. Trace Elements in Marine Sediment and Organisms in the Gulf of Thailand. *Int J Environ*
353 *Res Public Health* ~~[Internet]~~. 2018 ~~Apr 20~~ ~~[cited 2020 Aug 20]~~;15(4):810.
354 <https://doi.org/10.3390/ijerph15040810>
- 355 5. Joksimović D, Stanković S. Accumulation of trace metals in marine organisms of the southeastern Adriatic
356 coast, Montenegro. *J Serbian Chem Soc*. 2012 ~~Feb 2~~;77(1):105–17. <https://doi.org/10.2298/JSC110323159J>.
- 357 6. Chouvelon T, Strady E, Harmelin-Vivien M, Radakovitch O, Brach-Papa C, Crochet S, et al. Patterns of
358 trace metal bioaccumulation and trophic transfer in a phytoplankton-zooplankton-small pelagic fish marine
359 food web. *Mar Pollut Bull*. 2019 ~~Sep 1~~;146:1013–30. <https://doi.org/10.1016/j.marpolbul.2019.07.047>.
- 360 7. Shefer E, Silverman J, Herut B. Trace metal bioaccumulation in Israeli Mediterranean coastal marine

- 361 mollusks. *Quat Int.* 2015 ~~Dec-10~~;390:44–55. <https://doi.org/10.1016/j.quaint.2015.10.030>.
- 362 8. Galkus A, Joksas K, Stakeniene R, Lagunaviciene L. Heavy metal contamination of harbor bottom
363 sediments. *Polish J Environ Stud.* 2012;21(6):1583–94.
- 364 9. Kulkarni P, Chellam S, Ghurye G, Fraser MP. In Situ Generation of Hydrofluoric Acid during Microwave
365 Digestion of Atmospheric Particulate Matter Prior to Trace Element Analysis Using Inductively Coupled
366 Plasma Mass Spectrometry. *Environ Eng Sci* ~~[Internet]~~. 2003 ~~Jul-29~~ ~~[cited 2020 Aug-20]~~;20(6):517–31.
367 <https://www.liebertpub.com/doi/10.1089/109287503770736041>
- 368 10. Lo JM, Sakamoto H. Comparison of the Acid Combinations in Microwave-assisted Digestion of Marine
369 Sediments for Heavy Metal Analyses. *Anal Sci* ~~[Internet]~~. 2005 ~~[cited 2020 Aug-20]~~;21(10):1181–4. [http://](http://https://doi.org/10.2116/analsci.21.1181)
370 <https://doi.org/10.2116/analsci.21.1181>
- 371 11. Salomons W, Förstner U. Trace metal analysis on polluted sediments: Part II: Evaluation of environmental
372 impact. *Environ Technol Lett* ~~[Internet]~~. 1980 ~~[cited 2020 Aug-11]~~;1(11):506–17. Available from:
373 <https://www.tandfonline.com/doi/abs/10.1080/09593338009384007>
- 374 12. Chester R, Hughes MJ. A chemical technique for the separation of ferro-manganese minerals, carbonate
375 minerals and adsorbed trace elements from pelagic sediments. *Chem Geol.* 1967 ~~Jan-1~~;2(C):249–62.
376 [https://doi.org/10.1016/0009-2541\(67\)90025-3](https://doi.org/10.1016/0009-2541(67)90025-3).
- 377 13. Barber C. Major and trace element associations in limestones and dolomites. *Chem Geol.* 1974;14(4):273–
378 80. [https://doi.org/10.1016/0009-2541\(74\)90064-3](https://doi.org/10.1016/0009-2541(74)90064-3).
- 379 14. Usero J, Gamero M, Morillo J, Gracia I. Comparative study of three sequential extraction procedures for
380 metals in marine sediments. *Environ Int.* 1998 ~~May-1~~;24(4):487–96. [https://doi.org/10.1016/S0160-](https://doi.org/10.1016/S0160-4120(98)00028-2)
381 [4120\(98\)00028-2](https://doi.org/10.1016/S0160-4120(98)00028-2).
- 382 15. Gupta SK, Chen KY. Partitioning of trace metals in selective chemical fractions of nearshore sediments.
383 *Environ Lett.* 1975;10(2):129–58. <https://doi.org/10.1080/00139307509435816>.
- 384 16. Ure MU, Thomas R, Littlejohn D. Ammonium acetate extracts and their analysis for the speciation of metal
385 ions in soils and sediments. *Int J Environ Anal Chem* ~~[Internet]~~. 1993 ~~Jun-1~~ ~~[cited 2020 Aug-11]~~;51(1–

- 4):65–84. <https://www.tandfonline.com/doi/abs/10.1080/03067319308027612>.
17. Quevauviller P, Griepink B, Rauret G. Single and sequential extraction in sediments and soils. *Int J Environ Anal Chem.* 1993;51(1–4):231–5. <https://doi.org/10.1080/00139307509435816>.
18. Tessier A, Campbell PGC, Bisson M. Sequential Extraction Procedure for the Speciation of Particulate Trace Metals. *Anal Chem.* 1979;51(7):844–51. <https://doi.org/10.1021/ac50043a017>.
19. European Commission. The certification of the extractable contents (mass fractions) of Cd, Cr, Cu, Ni, Pb and Zn in freshwater sediment following a sequential extraction procedure - BCR 701 [Internet]. Community Research. 2001 [cited 2020 Aug 11]. p. 1–88. <https://op.europa.eu/en/publication-detail/-/publication/02c41803-8f6d-41f2-a73d-c47225ff8c6e>
20. López-Sánchez JF, Sahuquillo A, Fiedler HD, Rubio R, Rauret G, Muntau H, et al. CRM 601, A stable material for its extractable content of heavy metals. *Analyst* [Internet]. 1998 Jan 1 [cited 2020 Aug 20];123(8):1675–7. Available from: <https://pubs.rsc.org/en/content/articlehtml/1998/an/a802720j>
21. Kumkrong P, Mercier PHJ, Pihilligawa Gedara I, Mihai O, Tyo DD, Cindy J, et al. Determination of 27 metals in HISS-1, MESS-4 and PACS-3 marine sediment certified reference materials by the BCR sequential extraction. *Talanta* [Internet]. 2021 Jan 1 [cited 2020 Sep 17];221:121543. <https://doi.org/10.1016/j.talanta.2020.121543>
22. Sungur A, Soylak M, Ozcan H. Investigation of heavy metal mobility and availability by the BCR sequential extraction procedure: Relationship between soil properties and heavy metals availability. *Chem Speciat Bioavailab* [Internet]. 2014 [cited 2020 Aug 11];26(4):219–30. <https://www.tandfonline.com/doi/abs/10.3184/095422914X14147781158674>
23. Kubota R, Ohta A, Okai T. Speciation of 38 elements in eight GSJ geochemical sedimentary reference materials determined using a sequential extraction technique. *Geochem J* [Internet]. 2014 Mar 20 [cited 2020 Aug 10];48(2):165–88. <https://doi.org/10.2343/geochemj.2.0297> <http://jlc.jst.go.jp/DN/JST.JSTAGE/geochemj/2.0297?lang=en&from=CrossRef&type=abstract>

- 411 24. Long Y-Y, Hu L-F, Fang C-R, Wu Y-Y, Shen D-S. An evaluation of the modified BCR sequential
412 extraction procedure to assess the potential mobility of copper and zinc in MSW. *Microchem J.* 2009-Jan
413 4;91(1):1–5. <https://doi.org/10.1016/j.microc.2008.05.006>.
- 414 25. Craba L, Brunori C, Galletti M, Cremisini C, Morabito R. Comparison of three sequential extraction
415 procedures (original and modified 3 steps BCR procedure) applied to sediments of different origin. *Ann*
416 *Chim-[Internet]*. 2004 Jun 1 [cited 2020 Aug 11];94(56):409–19. Available from:
417 <http://doi.wiley.com/10.1002/adic.200490050>
- 418 26. Tokalioğlu Ş, Kartal Ş, Birol G. Comparison of three sequential extraction procedures for partitioning of
419 heavy metals in car park dusts. *J Environ Monit-[Internet]*. 2003 Jun 23 [cited 2020 Aug 11];5(3):468–76.
420 Available from: <https://pubs.rsc.org/en/content/articlehtml/2003/em/b300047h>
- 421 27. Berman S, Boyko V, Clancy V, Lam J, Maxwell P, McLaren J, et al. HISS-1 : Marine Sediment Reference
422 Material for Trace Metals and other Constituents [Internet]. 1997-[cited 2020 Aug 11]. [https://nrc-digital-](https://nrc-digital-repository.canada.ca/eng/view/object/?id=0f2ce620-b76b-4578-ba06-6aae6e632038)
423 [repository.canada.ca/eng/view/object/?id=0f2ce620-b76b-4578-ba06-6aae6e632038](https://nrc-digital-repository.canada.ca/eng/view/object/?id=0f2ce620-b76b-4578-ba06-6aae6e632038)
- 424 28. Willie S, Nadeau K, Gedara Pihillagawa I, Yang L, Clancy V, Grinberg P, et al. MESS-4 : Marine Sediment
425 Reference Material for Trace Metals and other Constituents [Internet]. 2014-[cited 2020 Aug 11]. Available
426 from: <https://nrc-digital-repository.canada.ca/eng/view/object/?id=8a3fd39a-c068-4ce0-820c-d08cf742a20a>
- 427 29. Willie S, Boyko V, Brophy C, Clancy V, Gedara Pihillagawa I, Grinberg P, et al. PACS-3 : Marine
428 Sediment Reference Material for Trace Metals and other Constituents [Internet]. 2013-[cited 2020 Aug 11].
429 Available from: [https://nrc-digital-repository.canada.ca/eng/view/object/?id=6616d284-f101-4b7a-b3c1-](https://nrc-digital-repository.canada.ca/eng/view/object/?id=6616d284-f101-4b7a-b3c1-a5280e4bd355)
430 [a5280e4bd355](https://nrc-digital-repository.canada.ca/eng/view/object/?id=6616d284-f101-4b7a-b3c1-a5280e4bd355)
- 431 30. Koepke A, Lafarge T, Possolo A, Toman B. NIST Consensus Builder—User’s Manual [Internet]. NIST,
432 editor. U.S. Department of Commerce; 2017.
433 https://consensus.nist.gov/app_direct/nicob/NISTConsensusBuilder-UserManual.pdf
- 434 31. Cox M, Harris P, Ellison S. Use of an “excess-variance” approach for the estimation of a key comparison
435 reference value, associated standard uncertainty and degrees of equivalence for CCQM key comparison
436 data. *CCQM/11-18-[Internet]*. 2011-[cited 2020 Aug 11];1–10.

- 437 <https://pdfs.semanticscholar.org/8f82/448ad8ae426308b1b631748d4e3ed69e891b.pdf>
- 438 32. Beldi G, Jakubowska N, Simoneau C. Report of the interlaboratory comparison organised by the European
439 Reference Laboratory for Food Contact Material - Formaldehyde in 3% acetic acid migration solution -
440 follow up exercise, EUR 26206 EN, [Internet]. 2013. [https://op.europa.eu/en/publication-detail/-](https://op.europa.eu/en/publication-detail/-/publication/6e9cdb24-f7e2-41c4-9edf-d21fe9daa9ae/language-en)
441 [/publication/6e9cdb24-f7e2-41c4-9edf-d21fe9daa9ae/language-en](https://op.europa.eu/en/publication-detail/-/publication/6e9cdb24-f7e2-41c4-9edf-d21fe9daa9ae/language-en)
- 442 33. Thompson M. The amazing Horwitz function [Internet]. amc technical brief. 2004 Jul [cited 2020 Aug 11].
443 Report No.: 17. https://www.rsc.org/images/horwitz-function-technical-brief-17_tcm18-214859.pdf
- 444 34. Horwitz W. Evaluation of Analytical Methods Used for Regulation of Foods and Drugs. Anal Chem
445 [Internet]. 1982 [cited 2020 Aug 20];54(1):67–76. [https://](https://pubs.acs.org/doi/abs/10.1021/ac00238a002)
446 pubs.acs.org/doi/abs/10.1021/ac00238a002 pubs.acs.org/sharingguidelines
- 447 35. Ellison SLR, Williams A, editors. Eurachem / CITAC Guide: Quantifying Uncertainty in Analytical
448 Measurement [Internet]. Third edit. Eurachem / CITAC Guide CG 4. Eurachem; 2012. 141 p.
449 https://www.eurachem.org/images/stories/Guides/pdf/QUAM2012_P1.pdf
- 450 36. Fleischer M. Recent estimates of the abundance of the elements in the earth's crust. Geol Survey Circ 285
451 [Internet]. 1953 [cited 2020 Aug 20];1–14. Available from: <https://pubs.er.usgs.gov/publication/cir285>
- 452 37. Filipkowska A, Kowalewska G, Pavoni B. Organotin compounds in surface sediments of the Southern
453 Baltic coastal zone: A study on the main factors for their accumulation and degradation. Environ Sci Pollut
454 Res. 2014;21(3):2077–87. <https://doi.org/10.1007/s11356-013-2115-x>.
- 455 38. Kelly WR, Murphy KE, Becker DA, Mann JL. Determination of Cr in certified reference material HISS-1,
456 marine sediment, by cold plasma isotope dilution ICP-MS and INAA: Comparison of microwave versus
457 closed (Carius) tube digestion. J Anal At Spectrom . 2003 Jan 28 [cited 2020 Aug 11];18(2):166–9.
458 <https://pubs.rsc.org/en/content/articlehtml/2003/ja/b210056h>
- 459

Figure

Tessier sequential extraction on 17 elements from three marine sediment certified reference materials (HISS-1, MESS-4 and PACS-3)

Paramee Kumkrong*, Ovi Mihai, Patrick H.J. Mercier, Indu Gedara Pihilligawa, Daniel Tyo, and Zoltan Mester

National Research Council Canada, 1200 Montreal Road Ottawa, Ontario Canada, K1A 0R6

*Corresponding author : Paramee.Kumkrong@nrc-cnrc.gc.ca




Figure 1 : The number of publications of Tessier and BCR sequential extraction procedure since 1979 to 2020
(Source: Scopus database, searched under keywords: “sequential extraction BCR” and “sequential extraction Tessier” (<https://www.scopus.com>, accessed on June 07, 2020))




Figure 2 : Tessier sequential extraction procedure flowchart




Figure 3 : Percent distribution (y-scale) of elements (x-scale) in Tessier extraction from A HISS-1, B MESS-4 and C PACS-3




Figure 4 : Element concentration and elements in HISS-1, MESS-4 and PACS-3 marine sediments, in comparison to Pelagic clay [3] and Earth crust [38]




Figure 5 : Standard uncertainty (% , at k=1) of A HISS-1, B MESS-4, and C PACS-3 from Step T # 1, T # 2, T # 3, T # 4 and residue in comparison with target uncertainty from Horwitz equation (red dot).




Figure 6 : Percent recovery of combined concentrations of four extraction steps and residue vs elements of A HISS-1, B MESS-4 and C PACS-3 with observed uncertainty (k=2)

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51
52
53
54
55
56
57
58
59
60

Figure captions

- Figure 1 : The number of publications of Tessier and BCR sequential extraction procedure since 1979 to 2020**
(Source: Scopus database, searched under keywords: “sequential extraction BCR” and “sequential extraction Tessier” (<https://www.scopus.com>, accessed on June 07, 2020))
- Figure 2 : Tessier sequential extraction procedure flowchart**
- Figure 3 : Percent distribution (y-scale) of elements (x-scale) by Tessier extraction from A HISS-1, B MESS-4 and C PACS-3**
- Figure 4 : Element concentration and elements in HISS-1, MESS-4 and PACS-3 marine sediments, in comparison to Pelagic clay [3] and Earth crust [38]**
- Figure 5 : Standard uncertainty (% , at k=1) of A HISS-1, B MESS-4, and C PACS-3 from Step T # 1, T # 2, T # 3, T # 4 and residue in comparison with target uncertainty from Horwitz equation (red dot).**
- Figure 6 : Percent recovery of combined concentrations of four extraction steps and residue vs elements of A HISS-1, B MESS-4 and C PACS-3 with observed uncertainty (k=2)**

Supplementary Material

Tessier sequential extraction on 17 elements from three marine sediment certified reference materials (HIS-1, MESS-4 and PACS-3)

Paramee Kumkrong*, Ovi Mihai, Patrick H.J. Mercier, Indu Gedara Pihilligawa, Daniel D. Tyo, and Zoltan Mester

National Research Council Canada, 1200 Montreal Road Ottawa, Ontario Canada, K1A 0R6

*Corresponding author : paramee.kumkrong@nrc-cnrc.gc.ca

Table S1 : Laboratory ware, chemicals and instruments

Laboratory ware	Lab A	Lab B
Extraction tube	A 50 mL polypropylene conical tube-Falcon with a screw cap from ThermoFisher Scientific Canada Pre-cleaned by soaking in 5% (v/v) nitric acid overnight, washing with double-distilled deionized water (later call water), and rinsed three times with water.	A 50 mL ultra-high performance centrifuge tube with a screw cap from VWR Canada. Pre-cleaned was conducted the same as Lab A
Chemicals	Lab A	Lab B
Deionized water (DIW)	Reverse osmosis connected to the ultrapure water type I by Millipore Synergy UV-R, Canada	Reverse osmosis connected to a Milli-Q Advantage A10, Canada.
Sodium acetate anhydrous (NaOAc Cas number 127-09-3)	ACS $\geq 99\%$, Sigma-Aldrich, Germany	$\geq 99\%$, Fisher Bio-reagent, Canada
100 % (v/v) acetic acid (AcOH, Cas number 64-19-7)	Suprapur, Merck, Canada	ACS Grade, VWR (BDH) Chemicals, Canada
Hydroxylammonium chloride (NH ₂ OH.HCl Cas number 5470-11-1)	ACS, Sigma-Aldrich, Canada	$\geq 96\%$ (w/w), Anachemia, Canada
30 % (w/w) hydrogen peroxide (H ₂ O ₂ , Cas number 7722-84-1)	Aristar Ultra, BDH, VWR Analytical, Canada	Optima Grade, ThermoFisher Chemicals Scientific, Canada
69 % (w/w) nitric acid (HNO ₃ , Cas number 7697-37-2)	Aristar Plus, BDH, VWR Analytical, Canada	J.T. Baker, Canada Distilled in house
ammonium acetate (NH ₄ OAc, Cas number 631-61-8)	$\geq 99\%$ (w/w) ACS, ThermoFisher Scientific, Canada.	$\geq 98\%$ (w/w) Sigma Aldrich, Canada.
48 % (w/v) hydrofluoric acid (HF, Cas number 7664-39-3)	Aristar Ultra, BDH, VWR Analytical, Canada.	Trace Metal Grade, ThermoFisher Chemical Scientific, Canada.
$\geq 99.5\%$ (w/w) boric acid (H ₃ BO ₃ , Cas number 10043-35-3)	Analar Grade, Chemique ACP Chemical Inc, Montreal Quebec Canada	Not use
37 % (w/w) hydrochloric acid (HCl, Cas number 7647-01-0)	Aristar Plus, BDH, VWR Analytical, Canada.	J.T. Baker, Canada, Distilled in house
Standard solution	- As, Co, Cr, Cu, Pb, Mn, Ni, Se, V, and Zn were ICP grade from SCP Science, Canada - Be, Cd, and Sb were ICP grade from ISO Spec, Delta, Canada - U was ICP grade from High Purity Standard, Canada. All standard solutions were traceable to SI unit via NIST.	CRMs from the NRC Canada for As, Cr, Mn, Ni, and Pb; and the rest of elements were based on CRMs from NIST.
Instruments	Lab A	Lab B

Analytical balance	Mettler Toledo AB304-S (Mettler Toledo GmbH, Laboratories & Weighing Technologies, Switzerland)	Mettler Toledo XPR305D5 (Mettler Toledo GmbH, Laboratories & Weighing Technologies, Switzerland)
Double distilled deionized water	Reverse osmosis connected to Millipore Synergy UV-R	Reverse osmosis connected to a Milli-Q Advantage A10
Shaker	End over end shaker Genie SI-1100 Roto-Shake Rotator (Scientific Industries Inc. USA)	Same as Lab A
Centrifugation*	Thermo Scientific Sorvall Legend XFR centrifuge at a speed of 3000 rpm or 2000xg (Thermo Fisher Scientific, Germany)	Thermo Scientific Sorvall Legend X1R Benchtop a speed of 3000 rpm or 1500xg (Thermo Fisher Scientific, Germany)
Hot block	Environmental express Model SC100 (Environmental express Inc., UA)	Same as Lab A
Microwave	Anton Paar Multiwave PRO with eight digestion tubes (Anton Paar GmbH, Austria)	Anton Paar Multiwave 3000 with 16 digestion tubes (Anton Paar GmbH, Austria)
Inductively coupled plasma-optical emission spectrometry (ICP-OES)**	Varian 720 (Agilent Technologies Inc, Australia)	- Agilent 5110 (SVDV) (Agilent Technologies Inc, USA), and - Varian 720 (Agilent Technologies Inc, Australia)
Inductively coupled plasma-mass spectrometry (ICP-MS)**	Agilent 7900 with collision reaction cell and MassHunter software (Agilent Technologies Inc, USA)	- Agilent 8800 triple Quadrupole ICP-MS (ICP-QQQ) with collision reaction cell and MassHunter software (Agilent Technologies Inc, USA), and - Perkin Elmer ELAN 9000 (PerkinElmer, Inc. USA)

*modified from the Tessier by using a speed of 3,000 rpm instead of 10,000 rpm (1,200 x g). The extraction solution after centrifugation was separated clearly from the sediment. Pasteur pipette or a narrow tip dropper was used to collect the extraction solution and to decant supernatant from the rinsing step. Care was taken not to disturb the sediments when the extraction solution level was close to the surface of sediments.

** Because of wide range concentrations of elements and matrix interferences, ICP-OES (axial and dual view) and various ICP-MS (single and triple quadrupole mass spectrometers) with a collision reaction cell were selected for comparison.

Table S2 : Instrumental conditions of elemental analysis

Condition	ICP-OES Varian 720	ICP-MS Agilent 7900	ICP-MS Agilent 8800
Radio frequency power (W)	1,200	1,550	1,560
Auxiliary argon flow rate (L/min)	0.7	1.0	0.9
Plasma gas flow (L/min)	12	15	15
Nebulizer flow (L/min)	0.7	1.04	1.1
Spray chamber	cyclonic	quartz Scott double pass	cyclonic
Sample pump rate (rpm)	12	18	15
Wavelength or mass	Nanometer (nm)	Isotopic mass	Isotopic mass
Arsenic (As)	193.696	75	75
Beryllium (Be)	*	9	9
Cadmium (Cd)	*	111	111
Chromium (Cr)	205.552	52	52 and 68 (O ₂)
Copper (Cu)	324.754	63	63
Iron (Fe)	259.940	-	-
Lithium (Li)	*	7	7
Manganese (Mn)	257.610	55	55
Nickel (Ni)	231.604	60	58 and 74 (O ₂) and 60 and 94 (NH ₃)
Lead (Pb)	220.353	206, 207, 208	208
Antimoni (Sb)	*	121, 123	121
Selenium (Se)	*	77, 78	82 (O ₂)
Tin (Sn)	*	118	118 and 135 (NH ₃)
Strontium (Sr)	421.552	-	88 and 104 (O ₂)
Uranium (U)	*	238	238 and 254 (O ₂)
Vanadium (V)	292.402	51	51 and 67 (O ₂)
Zinc (Zn)	213.856	66	64 and 80 (O ₂)
Auxiliary gas		No gas, H ₂ and He	No gas, O ₂ and NH ₃

*Analyzed by ICP-MS

Table S3 : Tessier sequential extraction procedure

Ext sol #	Reagent	Extraction Procedure
1	1 mol L ⁻¹ NaOAc, pH 8.2	-Extract 1 g of dried sample with 8 mL of Ext sol # 1, continuous end over end shaking at 30 ± 10 rpm for 1 h, centrifuge at 3000 rpm for 30 min, and collect the supernatant -Rinse the residue with 8 mL of DIW, shake for 15 min, centrifuge at 3000 rpm for 30 min, and discard the supernatant
2	1 mol L ⁻¹ NaOAc, pH 5	-Add 8 mL of Ext sol # 2, end over end shaking at 30 ± 10 rpm for 5 h, centrifuge at 3000 rpm for 30 min, and collect the supernatant - Rinse the residue with 8 mL of DI W, shake for 15 min, centrifuge at 3000 rpm for 30 min, and discard the supernatant
3	0.04 mol L ⁻¹ NH ₂ OH.HCl in 25 % (v/v) AcOH	-Add 20 mL of Ext sol # 3, place in 96 ± 3 °C hot block for 6 h with occasionally shaking, cool down, centrifuge at 3000 rpm for 30 min, and collect the supernatant - Rinse the residue with 8 mL of DIW, shake for 15 min, centrifuge at 3000 rpm for 30 min, and discard the supernatant
4	1) 0.02 mol L ⁻¹ HNO ₃ 2) 30 % H ₂ O ₂ at pH 2 3) 3.2 mol L ⁻¹ NH ₄ OAc in 20 % (v/v) HNO ₃	-Add 3 mL of 0.02 mol L ⁻¹ HNO ₃ and 5 mL of 30 % H ₂ O ₂ , stand for 15 min, place in 85 ± 2 °C hot block for 3 h with occasionally shaking and cool down -Add 3 mL of 30 % H ₂ O ₂ at pH 2, place in 85 ± 2 °C hot block for 3 h with occasionally shaking and cool down, -Add 5 mL of 3.2 mol L ⁻¹ NH ₄ OAc and 15 mL of DIW, end over end shaking at 30 ± 10 rpm for 30 min, centrifuge at 3000 rpm for 30 min, and collect the supernatant.

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47

Table S4 : Total concentrations and standard deviation, as mg kg⁻¹ dry weight, of elements in HISS-1, MESS-4 and PACS-3 as measured by two laboratories and the certified values where applicable

Elements	HISS-1 (mg kg ⁻¹)			MESS-4 (mg kg ⁻¹)			PACS-3 (mg kg ⁻¹)		
	Lab A	Lab B	Certified	Lab A	Lab B	Certified	Lab A	Lab B	Certified
	Found ± SD	Found ± SD	value ± U	Found ± SD	Found ± SD	value ± U	Found ± SD	Found ± SD	value ± U
As	0.863 ± 0.007	0.885 ± 0.031	0.801 ± 0.099	19.1 ± 2.0	23.6 ± 0.5	21.7 ± 2.8	24.6 ± 4.0	29.4 ± 1.0	30.3 ± 2.4
Be	0.110 ± 006	0.121 ± 0.010	0.129 ± 0.023	2.16 ± 0.07	1.93 ± 0.03	2.09 ± 0.28	1.03 ± 0.15	0.902 ± 0.102	1.06 ± 0.12
Cd	< 0.04	0.020 ± 0.006	0.024 ± 0.009	0.264 ± 0.042	0.316 ± 0.039	0.28 ± 0.04	2.49 ± 0.46	2.29 ± 0.16	2.23 ± 0.06
Cr	15.4 ± 2.2	15.9 ± 0.5	30 ± 6.8	86.7 ± 5.7	93.5 ± 3.5	94.3 ± 1.8	85.8 ± 6.3	89.5 ± 3.5	90.6 ± 4.0
Cu	1.97 ± 0.02	2.03 ± 0.25	2.29 ± 0.37	35.3 ± 5.8	30.1 ± 1.5	32.9 ± 1.8	292 ± 11	322 ± 7	326 ± 10
Fe	2,406 ± 203	2,363 ± 40	2,460 ± 90	36,067 ± 740	36,435 ± 1,770	37,900 ± 1,600	37,361 ± 1759	39,813 ± 199	41,060 ± 640
Li	2.80 ± 0.11	3.02 ± 0.05	2.83 ± 0.54	69.7 ± 6.0	76.8 ± 3.6	65.3 ± 6.8	36.9 ± 3.2	35.6 ± 0.8	31.9 ± 4.2
Mn	60.7 ± 4.0	63.7 ± 2.1	66.1 ± 4.2	289 ± 23	301 ± 15	298 ± 14	409 ± 28	424 ± 14	432 ± 16
Ni	2.55 ± 0.62	2.17 ± 0.21	2.16 ± 0.29	45.1 ± 4.2	41.5 ± 1.7	42.8 ± 1.6	41.7 ± 4.0	39.2 ± 2.0	39.5 ± 2.2
Pb	3.78 ± 0.88	3.00 ± 0.56	3.13 ± 0.40	19.8 ± 1.4	21.0 ± 0.9	21.5 ± 1.2	179.4 ± 16.1	189.4 ± 7.4	188.0 ± 7.4
Sb	0.092 ± 0.023	0.134 ± 0.009	0.13 *	1.13 ± 0.13	1.00 ± 0.16	1.07 ± 0.16	11.2 ± 1.4	12.0 ± 0.4	14.7 ± 2.2
Se	0.049 ± 0.009	0.056 ± 0.007	0.050 ± 0.007	0.963 ± 0.028	1.07 ± 0.04	1.5 *	0.955 ± 0.018	1.097 ± 0.186	Not reported
Sn	0.185 ± 0.006	0.132 ± 0.008	0.11 *	2.50 ± 0.11	2.50 ± 0.12	2.35 ± 0.12	26.4 ± 1.1	20.5 ± 1.7	22.0 ± 0.2
Sr	92.3 ± 3.3	100.6 ± 1.4	96.9 ± 11.2	115 ± 3	133 ± 4	132 ± 8	230 ± 3	268 ± 6	267 ± 10
U	0.352 ± 0.081	0.301 ± 0.019	0.26 *	3.20 ± 0.42	3.70 ± 0.22	3.4 ± 0.4	2.55 ± 0.32	2.59 ± 0.08	2.6 *
V	5.99 ± 0.22	7.22 ± 0.31	6.8 ± 0.78	197 ± 9	211 ± 5	216 ± 8	127 ± 4	129 ± 4	129 ± 8
Zn	6.87 ± 1.82	4.84 ± 0.45	4.94 ± 0.79	134 ± 17	142 ± 6	147 ± 6	352 ± 31	374 ± 6	376 ± 12

* an indicative value




Figure S1: Percent recovery (y-scale) of total element (x-scale) determination by two laboratories in comparison with certified values




Figure S2 : Distribution histogram of observed standard uncertainties of elements in HISS-1, MESS-4 and PACS-3 (bar chart) and target uncertainty (line plot) by Horwitz equation

Table S5 : Combined concentrations (mg kg^{-1} dry weight) of elements from four extraction steps and residue and associated uncertainties (U, $k=2$) in comparison to certified values with associated uncertainties ($k=2$) of HISS-1 and percent recovery

Metals	Combined value	U	Certified value	Certified U	Recovery (%)
Li	3.60	2.64	2.83	0.54	127
Be	0.125	0.059	0.129	0.023	97
V	8.06	2.12	6.80	0.78	119
Cr	18.5	0.6	30.0	6.8	62
Mn	67.0	7.1	66.1	4.2	101
Fe	2,361	13	2,460	90	96
Ni	2.91	0.33	2.16	0.29	135
Cu	2.479	0.324	2.290	0.370	108
Zn	5.07	0.73	4.94	0.79	103
As	0.902	0.407	0.801	0.099	113
Sr	82.9	18.2	96.9	11.2	86
Sn	0.192	0.114	0.142	0.020	135
Sb	0.116	0.005	0.119	0.049	98
Pb	3.48	0.29	3.13	0.40	111
U	0.228	0.043	0.304	0.166	75
Total concentration of U used average concentration of two laboratories					

Table S6 : Combined concentrations (mg kg⁻¹ dry weight) of elements from four extraction steps and residue and associated uncertainties (U, k=2) in comparison to certified values with associated uncertainties (k=2) of MESS-4 and percent recovery

Metals	Combined value	U	Certified value	Certified U	Recovery (%)
Li	58.7	7.9	65.3	6.8	90
Be	2.30	0.22	2.09	0.28	110
V	222	25	216	8	103
Cr	94.5	10.5	94.3	1.8	100
Mn	273	20	298	14	91
Fe	34818	6761	37900	1600	92
Ni	45.8	5.8	42.8	1.6	107
Cu	34.4	2.9	32.9	1.8	105
Zn	157	13	148	6	106
As	21.2	3.7	21.7	2.8	98
Se*	1.074	0.274	1.179	0.101	91
Sr	124	14	132	8	94
Cd	0.241	0.103	0.280	0.040	86
Sn	3.23	0.48	2.35	0.12	137
Sb	1.172	0.182	1.070	0.160	110
Pb	21.3	3.3	21.5	1.2	99
U	3.09	0.95	3.40	0.40	91
*Total concentration of Se used an average concentration of two laboratories					

Table S7 : Combined concentrations (mg kg⁻¹ dry weight) of elements from four extraction steps and residue and associated uncertainties (U, k=2) in comparison to certified values with associated uncertainties (k=2) of PACS-3 and percent recovery

Metals	Combined value	U	Certified value	Certified U	Recovery (%)
Li	30.2	5.0	31.9	4.2	95
Be	0.99	0.06	1.06	0.12	93
V	139	20	129	8	108
Cr	99.6	9.0	90.6	4.0	110
Mn	396	39	432	16	92
Fe	37086	870	41060	640	90
Ni	41.4	5.0	39.5	2.2	105
Cu	318	23	326	10	98
Zn	430	44	376	12	114
As	30.8	3.9	30.3	2.4	102
Se*	1.18	0.25	1.03	0.26	115
Sr	232	103	267	10	87
Cd	2.27	0.50	2.23	0.16	102
Sn	26.3	7.9	22.0	2.20	120
Sb	13.8	1.4	14.7	2.200	94
Pb	216	30	188	7.4	115
U*	2.64	0.45	2.58	0.66	102
*Total concentrations of Se and U used an average concentration of two laboratories					

Graphical Abstract

