



NRC Publications Archive Archives des publications du CNRC

Annotating and quantifying pri-miRNA transcripts using RNA-Seq data of wild type and serrate-1 globular stage embryos of Arabidopsis thaliana

Lepe-soltero, Daniel; Armenta-medina, Alma; Xiang, Daoquan; Datla, Raju; Gillmor, C. Stewart; Abreu-goodger, Cei

This publication could be one of several versions: author's original, accepted manuscript or the publisher's version. / La version de cette publication peut être l'une des suivantes : la version prépublication de l'auteur, la version acceptée du manuscrit ou la version de l'éditeur.

For the publisher's version, please access the DOI link below. / Pour consulter la version de l'éditeur, utilisez le lien DOI ci-dessous.

Publisher's version / Version de l'éditeur:

<https://doi.org/10.1016/j.dib.2017.10.019>

Data in Brief, 15, pp. 642-647, 2017-10-12

NRC Publications Record / Notice d'Archives des publications de CNRC:

<https://nrc-publications.canada.ca/eng/view/object/?id=6731d37e-5dca-4194-944e-dc8bc9ae247e>

<https://publications-cnrc.canada.ca/fra/voir/objet/?id=6731d37e-5dca-4194-944e-dc8bc9ae247e>

Access and use of this website and the material on it are subject to the Terms and Conditions set forth at

<https://nrc-publications.canada.ca/eng/copyright>

READ THESE TERMS AND CONDITIONS CAREFULLY BEFORE USING THIS WEBSITE.

L'accès à ce site Web et l'utilisation de son contenu sont assujettis aux conditions présentées dans le site

<https://publications-cnrc.canada.ca/fra/droits>

LISEZ CES CONDITIONS ATTENTIVEMENT AVANT D'UTILISER CE SITE WEB.

Questions? Contact the NRC Publications Archive team at

PublicationsArchive-ArchivesPublications@nrc-cnrc.gc.ca. If you wish to email the authors directly, please see the first page of the publication for their contact information.

Vous avez des questions? Nous pouvons vous aider. Pour communiquer directement avec un auteur, consultez la première page de la revue dans laquelle son article a été publié afin de trouver ses coordonnées. Si vous n'arrivez pas à les repérer, communiquez avec nous à PublicationsArchive-ArchivesPublications@nrc-cnrc.gc.ca.





ELSEVIER

Contents lists available at ScienceDirect

Data in Brief

journal homepage: www.elsevier.com/locate/dib

Data Article

Annotating and quantifying pri-miRNA transcripts using RNA-Seq data of wild type and *serrate-1* globular stage embryos of *Arabidopsis thaliana*

Daniel Lepe-Soltero^a, Alma Armenta-Medina^a,
 Daoquan Xiang^b, Raju Datla^b, C. Stewart Gillmor^a,
 Cei Abreu-Goodger^{a,*}

^a Laboratorio Nacional de Genómica para la Biodiversidad (Langebio), Unidad de Genómica Avanzada, Centro de Investigación y de Estudios Avanzados (CINVESTAV), Irapuato, Guanajuato, Mexico

^b Plant Biotechnology Institute, National Research Council, Saskatoon, Saskatchewan, Canada

ARTICLE INFO

Article history:

Received 13 September 2017

Received in revised form

5 October 2017

Accepted 6 October 2017

Available online 12 October 2017

ABSTRACT

The genome annotation for the model plant *Arabidopsis thaliana* does not include the primary transcripts from which MIRNAs are processed. Here we present and analyze the raw mRNA sequencing data from wild type and *serrate-1* globular stage embryos of *A. thaliana*, ecotype Columbia. Because *SERRATE* is required for pri-miRNA processing, these precursors accumulate in *serrate-1* mutants, facilitating their detection using standard RNA-Seq protocols. We first use the mapping of the RNA-Seq reads to the reference genome to annotate the potential primary transcripts of MIRNAs expressed in the embryo. We then quantify these pri-miRNAs in wild type and *serrate-1* mutants. Finally, we use differential expression analysis to determine which are up-regulated in *serrate-1* compared to wild type, to select the best candidates for *bona fide* pri-miRNAs expressed in the globular stage embryos. In addition, we analyze a previously published RNA-Seq dataset of wild type and *dicer-like 1* mutant embryos at

DOI of original article: <http://dx.doi.org/10.1016/j.ydbio.2017.09.009>

* Corresponding author.

E-mail address: cei.abreu@cinvestav.mx (C. Abreu-Goodger).

<http://dx.doi.org/10.1016/j.dib.2017.10.019>

2352-3409/© 2017 The Authors. Published by Elsevier Inc. This is an open access article under the CC BY license (<http://creativecommons.org/licenses/by/4.0/>).

the globular stage [1]. Our data are interpreted and discussed in a separate article [2].

© 2017 The Authors. Published by Elsevier Inc. This is an open access article under the CC BY license (<http://creativecommons.org/licenses/by/4.0/>).

Specifications Table

Subject area	Biology
More specific subject area	Plant biology
Type of data	Tables
How data was acquired	RNA-Seq from an Illumina HiSeq 2000 and previously published <i>dcl1-5</i> raw data from the Gene Expression Omnibus (GEO) accession GSE25404 [1]
Data format	Raw and analyzed data
Experimental factors	Total RNA was extracted from pools of ~80 globular stage (32–64 cell) embryos isolated at 72 hours after pollination
Experimental features	<i>Arabidopsis thaliana</i> wt and <i>se-1</i> embryos at the globular stage, with two biological replicates
Data source location	Not applicable
Data accessibility	Data is available as Supplementary file 1 and at NCBI GEO accession GSE100450

Value of the data

- This is the first study to directly identify *MIRNA* genes expressed in early embryos of plants.
- We provide an annotation file with 318 *MIRNA* gene models, including 77 predicted from the RNA-Seq data, that is useful for others interested in *MIRNA* gene regulation in *Arabidopsis*.
- Our high-quality globular stage transcriptomes of wild type and *serrate-1* embryos will be valuable for other studies of gene regulation in early embryogenesis.

1. Data

We generated RNA-Seq data for globular stage *Arabidopsis thaliana* embryos from two genotypes: *serrate-1* (*se-1*) mutants, and wild type (wt), both in the Columbia ecotype. We then inferred *MIRNA* primary transcripts expressed at the globular stage by aligning RNA-Seq reads to the *Arabidopsis* genome, assembling and manually curating gene models, and analyzing differential expression. As an independent profile of *MIRNA* transcripts in embryos, we analyzed a previously published RNA sequencing experiment using wt Columbia and *dicer-like 1* (*dcl1-5*) embryos [1]. We provide the raw data, predicted pri-miRNA gene models, quantification of all genes in both experiments, and differential expression results.

2. Experimental design, materials and methods

Two biological replicates of wt and *se-1* of about 80 embryos each at the 32–64 cell (early to mid-globular) stage were obtained, the RNA isolated, amplified and sequenced as described previously [2]. Illumina HiSeq 2000 sequencing yielded 101 nt paired-end reads with over 20 million reads per library (Table 1). Raw data files are available through the NCBI Gene Expression Omnibus (GEO, accession [GSE100450](#)). Raw data from the *dcl1-5* RNA-Seq experiment was downloaded from GEO

(accession GSE25404), consisting of a single replicate of *dcl1-5* mutant and Col-0 (wt) early globular embryos with 36 nt reads obtained on an Illumina Genome Analyzer II [1].

The paired-end reads from the wt and *se-1* libraries were mapped using HISAT2 [3] with default settings except for an intronic length suited for Arabidopsis (`-max-intronlen 900`). The default intronic length (`-max-intronlen 500000`), tuned for mammalian genomes, resulted in many reads falsely mapping across several genes. The 36 nt reads from the wt and *dcl1-5* libraries were cleaned using cutadapt [4] v1.13 to remove adaptors, polyA and polyT sequences from the 3' and 5' ends (respectively), and low-quality bases and flanking Ns were trimmed from individual reads. Finally, any read shorter than 18 nt or with more than three internal Ns was discarded. The full parameters for cutadapt were: `-q 6 -a ATCTCGTATGCCNNNNNNNNNNNNNNNNNNNNNNNNNNNNNN -a "A{36}" -g NNNNNNNNNNNNNNNNNNNNNNNNNNNNNAGTCC-GACGATC -g "T{36}" -trim-n -max-n 3 -m 18`.

The resulting cleaned reads were mapped using the Bowtie [5] short read aligner allowing up to 2 mismatches within a 25 nt seed sequence, and only uniquely mapped reads were retained (`-1 25 -n 2 -m 1`) (Table 2). Both transcriptomes were mapped to the reference TAIR10 assembly of *Arabidopsis thaliana*, downloaded from The Arabidopsis Information Resource (ftp://ftp.arabidopsis.org/home/tair/Sequences/whole_chromosomes/).

The Araport11 reference annotation from the Arabidopsis Information Portal [6] contains the coordinates of pre-miRNA hairpins, but does not contain information regarding the pri-miRNA transcripts. In order to predict genome wide pri-miRNA transcript coordinates, Cufflinks [7] was used to assemble and merge putative pri-miRNA transcripts from the *se-1* and wt Col RNA-Seq read alignments. Cufflinks was first run for each library (`-overlap-radius 1 -library-type fr-unstranded`). The predictions were then merged (`cuffmerge -s`) using the TAIR10 genome assembly as reference. Out of the 325 miRNAs in the Araport11 annotation, 77 overlapped with a predicted Cufflinks gene model. All these predictions were manually verified, with only 4 of them requiring individual adjustments to better reflect the disposition of the reads from the RNA-Seq libraries, and to resolve overlapping conflicts. The main limitation of this approach is that pri-miRNAs can only be assembled if they are expressed in the sampled conditions (in this case, for early globular embryos).

Overlaps with, or even proximity to protein-coding genes can make it difficult to establish the appropriate gene model of a pri-miRNA. Due to this, the pri-miRNA predictions were divided into four groups: intergenic (G1), between 1–400 bp away from a protein-coding gene (G2), overlapping with a protein-coding gene (G3, divided into G3A if the overlap includes the pre-miRNA or G3B otherwise) and overlapping with a non-coding gene (G4); see Table 3. The 77 Cufflinks predictions were distributed amongst the pri-miRNA groups as follows: 28 in G1, 7 in G2, 0 in G3A, 37 in G3B and 5 in G4. There were 56 pre-miRNAs that overlapped with 54 protein-coding genes (G3A group). These pre-miRNAs were assigned to a pri-miRNA gene model identical to the overlapping protein-coding gene.

Table 1
RNA sequences obtained by Illumina sequencing.

Sample	Raw reads	Quantified reads
WT1	20,751,215	18,894,161 (91.05%)
WT2	21,623,459	19,362,653 (89.54%)
<i>se1</i>	20,708,775	18,580,516 (89.72%)
<i>se2</i>	22,984,350	20,882,004 (90.85%)

Table 2
RNA sequences from GSE25404.

Sample	Raw reads	Quantified reads
WT (Col-0) (SRR074122)	21,413,867	5,153,766 (24.06%)
<i>dcl1-5</i> (SRR074123)	22,578,840	4,047,267 (17.92%)

Table 3

Groups of miRNAs according to their position relative to other annotation in the Arabidopsis genome.

Group	Description	Number of pre-miRNAs in Araport11	Predicted pri-miRNA gene models	Gene models from Araport11	Gene models with two pre-miRNAs	Pri-miRNAs gene models in final annotation
G1	Intergenic region	222	28	166	2	194
G2	Closest protein-coding gene is 1–400 bp away	39	7	19	0	26
G3A	Overlap of pre-miRNA with protein-coding gene	56	0	54 ^a	3	54
G3B	Overlap of pri-miRNA with protein-coding gene	0	37	0	1	37
G4	Overlaps with non-coding gene	8	5	2	1	7
Total		325	77	241	7	318

^a The overlapping protein gene model was used instead of the pre-miRNA coordinates.**Table 4**Common pri-miRNAs detected with at least 4 reads in the *se-1* and *dcl1-5* experiments.

	<i>dcl1-5</i> (121 miRNAs)
<i>se-1</i> (133 miRNAs)	117
<i>se-1</i> (FDR < 0.05) (100 miRNAs)	90

In all other cases where no Cufflinks prediction was available, the pri-miRNA was kept the same as the pre-miRNA annotation from Araport11.

A final annotation file with the newly predicted pri-miRNA gene models, in addition to all the gene models from Araport11, considering a total of 318 pri-miRNA genes and 27,562 protein-coding genes, was employed for the quantification of all the RNA-Seq libraries and is available as [Supplementary file 1](#). Quantification of reads using this annotation file was performed in R with the function *featureCounts* from the *Rsubread* package [8]. Multi-mapping reads were counted (`countMultiMappingReads=TRUE`) and only primary alignments were allowed (`primaryOnly=TRUE`). Additionally, reads were assigned to the feature with the largest number of overlapping bases (`largestOverlap=TRUE`) and a minimum mapping quality score of 10 was required (`minMQS=10`) for a read to be counted.

Finally, the *edgeR* [9] package was used to perform the differential expression analysis of both the *se* and *dcl1* experiments, using the raw counts with no prior filtering. A tagwise dispersion was calculated for *se*, but since no replicates are available for the *dcl1* experiment, the Biological Coefficient of Variation was fixed to 0.4, as recommended by the *edgeR* manual. To test for differential expression, quasi-likelihood F-tests and likelihood ratio tests were performed for the *se* and *dcl1* experiments, respectively. In total, 6951 genes were upregulated in the *se-1* mutant and 7138 were downregulated with an FDR < 0.05, and 125 genes were upregulated in *dcl1-5* and 138 downregulated with an FDR < 0.05 (Table S3 from reference [2]).

Of the 318 annotated pri-miRNAs, 100 were deemed differentially expressed (FDR < 0.05) in *se-1*. Of those, 73 were upregulated (G1: 22, G2: 6, G3A: 17, G3B: 25, G4: 3) and 27 were downregulated (G1: 5, G2: 1, G3A: 14, G3B: 7, G4: 0). Because *SERRATE* participates during pri-miRNA processing, they should be up-regulated in *se-1*. To evaluate which of the pri-miRNA groups behaved as expected in the *se-1* mutant, a one-sided Wilcoxon rank sum test was done with the *wilcox.test* function from the *stats* package in R. The resulting p-values were G1: 0.0016, G2: 0.017, G3A: 0.14, G3B: 0.000017, and G4: 0.04, suggesting that in most cases the gene-models do reflect the properties of a pri-miRNA

Table 5
Behavior of differentially expressed pri-miRNAs in *se-1* and *dcl1-5*.

117 common pri-miRNAs (no FDR filter)		<i>dcl1-5</i>	
		Upregulated	Downregulated
<i>se-1</i>	Upregulated	47	39
	Downregulated	11	20

Table 6
Behavior of differentially expressed pri-miRNAs in *se-1* (FDR < 0.05) and *dcl1-5* (no FDR filter).

90 common pri-miRNAs (FDR < 0.05 for <i>se-1</i>)		<i>dcl1-5</i>	
		Upregulated	Downregulated
<i>se-1</i>	Upregulated	37	31
	Downregulated	7	15

(except for the G3A category, where the gene models are taken from protein-coding genes). The average values of expression and \log_2FC for all genes, including pri-miRNAs and the median \log_2FC for each of the miRNA groups, is plotted in Fig. S2 from Ref. [2].

In *se-1*, 133 miRNAs have at least 4 accumulated reads from the libraries. Of these, 100 are differentially expressed (FDR < 0.05). For *dcl1-5*, 121 miRNAs have at least 4 accumulated reads from the libraries. A summary of the miRNAs that were detected in both the *se-1* and *dcl1-5* experiments is detailed in Table 4 and a summary of their differential expression behavior is shown in Tables 5 and 6.

Acknowledgements

This study was financed by CINVESTAV institutional funds provided to CSG, and by a HFSP Young Investigator Award (RGY0069) to CA-G.

Transparency document. Supporting information

Transparency data associated with this article can be found in the online version at <http://dx.doi.org/10.1016/j.dib.2017.10.019>.

Appendix A. Supporting information

Supplementary data associated with this article can be found in the online version at <http://dx.doi.org/10.1016/j.dib.2017.10.019>.

References

- [1] M.D. Nodine, D. Bartel, MicroRNAs prevent precocious gene expression and enable pattern formation during plant embryogenesis, *Genes Dev.* 24 (2010) 2678–2692.
- [2] A. Armenta-Medina, D. Lepe-Soltero, D. Xiang, R. Datla, C. Abreu-Goodger, C.S. Gillmor, Arabidopsis thaliana miRNAs promote embryo pattern formation beginning in the zygote, *Dev. Biol.* (2017), <http://dx.doi.org/10.1016/j.ydbio.2017.09.009>.
- [3] D. Kim, B. Langmead, S.L. Salzberg, HISAT: a fast spliced aligner with low memory requirements, *Nat. Methods* 12 (2015) 357–360. <http://dx.doi.org/10.1038/nmeth.3317>.

- [4] M. Martin, Cutadapt removes adapter sequences from high-throughput sequencing reads, *EMBnet.J.* 17 (2011) 10. <http://dx.doi.org/10.14806/ej.17.1.200>.
- [5] B. Langmead, C. Trapnell, M. Pop, S.L. Salzberg, Ultrafast and memory-efficient alignment of short DNA sequences to the human genome, *Genome Biol.* 10 (2009) R25. <http://dx.doi.org/10.1186/gb-2009-10-3-r25>.
- [6] C.-Y. Cheng, V. Krishnakumar, A.P. Chan, F. Thibaud-Nissen, S. Schobel, C.D. Town, Araport11: a complete reannotation of the *Arabidopsis thaliana* reference genome, *Plant J.* 89 (2017) 789–804. <http://dx.doi.org/10.1111/tpj.13415>.
- [7] C. Trapnell, B.A. Williams, G. Pertea, A. Mortazavi, G. Kwan, M.J. van Baren, S.L. Salzberg, B.J. Wold, L. Pachter, Transcript assembly and quantification by RNA-Seq reveals unannotated transcripts and isoform switching during cell differentiation, *Nat. Biotechnol.* 28 (2010) 511–515.
- [8] Y. Liao, G.K. Smyth, W. Shi, FeatureCounts: an efficient general purpose program for assigning sequence reads to genomic features, *Bioinformatics* 30 (2014) 923–930. <http://dx.doi.org/10.1093/bioinformatics/btt656>.
- [9] M.D. Robinson, D.J. McCarthy, G.K. Smyth, edgeR: a bioconductor package for differential expression analysis of digital gene expression data, *Bioinformatics* 26 (2010) 139–140.